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**PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

### Bivalve cervical collar with open trachea

#### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

#### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

#### SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields. Do not apply in direct contact with open wounds.

#### WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

#### SELECTION/SIZE

Code	REF. 1530/T		
Size	S	M	L
Neck circ. cm	30/35	35/40	40/45
Height neck* cm	2 <sup>1/4"</sup> = 6 cm · 3 <sup>1/4"</sup> = 8 cm · 4 <sup>1/4"</sup> = 11 cm · 5 <sup>1/4"</sup> = 13 cm		
Colour	beige		

\* vertical distance between the upper side of the shoulder and the chin

#### MAINTENANCE

<input checked="" type="checkbox"/> Do not bleach	<input checked="" type="checkbox"/> No chemical cleaning
<input checked="" type="checkbox"/> Do not iron	<input checked="" type="checkbox"/> Do not tumble-dry
<input checked="" type="checkbox"/> Washing instructions: hand wash at max. 30° C with neutral soap (a sponge is recommended); do not dry in the vicinity of heat sources.	

Do not dispose of the product or any of its components into the environment.

#### INDICATIONS

- Severe to serious distortions of cervical spine (whiplash)
- Multiple disc diseases of the cervical spine
- Cervical trauma
- Pathological failure due to cervical metastases
- Use after cervical spine surgery

#### CONTRAINDICATIONS

Currently no known

#### CHARACTERISTICS AND MATERIALS

- Structure made of wrapping bivalve with Back Velcro® fastening system
- Structure made of non-citotoxic polymeric foam in compliance with UNI EN ISO 10993-5
- Rigid plates for back, front and central support
- Radiolucent and waterproof
- Shaped for a perfect adhesion to chin, neck and shoulders
- Open trachea and holes for ventilation
- Prepared to be connected with chest stabilizer or hyperextensio

#### PUTTING ON THE APPLIANCE

- Position the front shell in such a way that the arrow points upwards. The chin should be in the cut-out area (fig. A)
- Position the rear shell in such a way that the arrow points upwards. Place the side loops over one another. In doing so, those on the rear shell should be inside and those on the front shell should be outside (fig. B)
- Pull both loops evenly at both sides and secure well using the Velcro® closures (fig. C)
- If necessary, the rigid support plates in the middle may be heated to be bent into the desired form. Any protruding foam body elements may be cut to size using scissors.

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