 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Dorsal-sternum stabilizer for Cervistable collar

DECLARATION OF CONFORMITY

As the distributor, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.1521
Size	universal
Colour	beige

MAINTENANCE

Washing instructions: wipe the internal and external parts with a soft sponge dipped in lukewarm water and mild soap. Leave to dry away from direct heat.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Cervical-dorsal fractures
- Cervical spondylolisthesis
- Multiple discopathy
- Cervical and breastbone fractures due to osteoporosis or metastasis
- Use after spine surgery and as an immobilizer in the event of first aid

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- PVC struts
- PLASTAZOTE® cover
- Connections with fixing buttons to Cervistable collar
- Velcro® fastening system - lower part
- Radiolucent

APPLICATION

Proper application of the Stabilizer is important for immobilization and patient comfort. The two-piece design consists of a FRONT piece and BACK piece which are packaged as a set.

Always maintain the patient's head in neutral alignment. These directions are guidelines only; follow a qualified healthcare provider's instructions.

1. Insert the plastic reinforcement of the Collar FRONT piece into the top of the Stabilizer FRONT piece. Slide the pieces together until the holes are in alignment. Repeat for the Collar BACK piece and Stabilizer BACK piece.
2. After the Stabilizer and Collar holes are in alignment, insert the locking plugs. Detach the nylon fastener after the plugs are securely locked into position and the two pieces cannot be separated.
3. After properly applying the Collar, insert the Velcro straps located on the Stabilizer BACK piece through the buckles located on the Stabilizer FRONT piece.
4. Pull both straps simultaneously and with equal pressure toward the Stabilizer BACK piece. Verify that all components are centered on the patient and the patient is in neutral alignment. Fasten strap securely and remove excess material, as necessary.

REMOVAL:

Before removing the Stabilizer and Collar, mark the position of the Velcro strap location to ensure re-application maintains the proper fit.

With the Stabilizer and Collar secured on the patient, loosen and unbuckle the Stabilizer strap on the RIGHT side and then on the LEFT side. To remove the Collar, unfasten the Velcro strap on the RIGHT side and then the LEFT side. Remove Stabilizer and Collar from the patient.

Plugs can be removed only when the Stabilizer is no longer necessary and only the Collar is required as instructed by a qualified healthcare professional.

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ORTHOSERVICE DEUTSCHLAND GmbH
OSUK Ltd
RO+TEN S.r.l.