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SINGLE PATIENT MULTIPLE USE

LOW PROFILE

 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

## Back orthopedic brace

### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

### SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields. Do not apply in direct contact with open wounds.



### WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

### SELECTION/SIZE

Item	REF.91111
Size	universal
Height cm	30
Colour:	cream

### MAINTENANCE

-  Do not bleach  No chemical cleaning
-  Do not iron  Do not tumble-dry
-  Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

### INDICATIONS

- Useful for the correction of a curved back, kyphosis and scapular winging.
- A useful aid for back pain caused by poor posture.

### CONTRAINDICATIONS

Currently no known

### CHARACTERISTICS AND MATERIALS

- The support and straps are made of PE foam in a bi-layer of velveteen and technical materia
- A back pulley system provides better traction when fastening.
- The straps can be adjusted and shortened using the Y-Velcro® straps:
  - The overall length can be adjusted;
  - The positioned of the shoulder straps can be adjusted.
- Fastens with a snap-fit buckle.
- Easy to put on, just like a backpack.

### PUTTING ON THE APPLIANCE

#### WEARING THE BRACE FOR THE FIRST TIME

- To put on the brace slide the patient’s arms through the shoulder straps (fig A).
- Adjust the position of the shoulder straps using the Y-tab on the top end (fig B):
  - Position the shoulder straps so that the padding protects the armpit.
  - If they are too long, remove the Y-Velcro® strap, cut it to the correct length and reposition the Y-Velcro® strap.
- Fasten the snap-fit buckle at the front (fig C).
- Adjust the length of the lower straps:
  - Stretch the straps to fit the patient (fig D).
  - If they are too long, remove the Y-Velcro® strap, cut it to the correct length and reposition the Y-Velcro® strap (fig H + fig I).
- Place the brace back on the patient and adjust the fit if necessary.

### WEARING THE BRACE ONCE FITTED

- Unfasten the snap-fit buckle.
- To put on the brace slide the patient’s arms through the shoulder straps (fig E).
- Pull the straps at the sides (fig F) and fasten them at the front with the snap-fit buckle (fig G).

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