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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

REF.55113 Post-operative thoracic band support REF.55114 Post-operative lumbar band

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/ or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.


WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.55113 • REF.55114					
Size	XS	S	M	L	XL	XXL
Chest-waist circumference cm	<100	100/110	110/130	130/140	140/150	150/160
Height cm REF.55113	16	16	16	16	16	16
Height cm REF.55114	24	24	24	24	24	24
Colour	white					

MAINTENANCE

 Do not bleach  No chemical cleaning
 Do not iron  Do not tumble-dry
 Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Post-operative hospital care
- Only REF.55113**
- Thoracic contusions
- Only REF.55114**
- Umbilical hernias
- Post-natal containment support
- Eventration
- Abdominal hernias

CONTRAINDICATIONS

There are no known contraindications at this time.

FEATURES AND MATERIALS

- Breathable
- Front Velcro® fastening

Only REF.55113

- The two elasticated bands allow the belt to adapt and fit well on the body

Only REF.55114

- The three elastic bands facilitate a good fit which easily adapts to the body

PUTTING ON THE APPLIANCE

Only REF.55113

- Wrap the support around the thorax and fasten with the Velcro® straps provided

Only REF.55114

- Wrap the support around the abdomen and fasten with the Velcro® straps provided.

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