MOLIMO VAS DA PAŽLJIVO PROČITATE OVE UPUTE I SAČUVATE IH
Elastična ortaza s rigidnim i termoformabilnim ispunjenjem i sustavom remenica za stabilizaciju lumbalne kralježnice i držanja

IZJAVA O SUKLADNOSTI

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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Elastic orthosis with rigid and thermoformable pads to stabilize the lumbar rachis and relieve pressure

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V. Sensitive® is a registered trademark of EUROJERSEY

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not too tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the device, to avoid the risk of irritation. Please refer to the instructions for use of the device. In the event of intolerance or allergies to one or more components in the product. We recommend not to wear the device in the proximity of fire flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist or an orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices in hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF. 1090 - REF. 1091*
Size	XS S M L XL XXL
Pelvic circumf. cm	70/80 80/90 90/100 100/110 110/120 120/130
Front height cm	19 19 19 19 19 19
Front height cm*	23 23 23 23 23 23
Back height cm	28 28 28 28 28 28
Back height cm*	33 33 33 33 33 33
Colour	C-cream (1090,1091) / S-black (1090,1091)

OPTION 95T012/L optional rigid stays cm 22 (REF.1090) - cm 26 (REF.1091)

MAINTENANCE
☒ Do not bleach ☒ No chemical cleaning
☒ Do not iron ☒ Do not tumble-dry
☒ Washing instructions:

ORTHOSES

- ensure that the Velcro® fasteners are securely closed before washing the orthosis
- wash the orthosis separately
- wash by hand or in the washing machine with mild detergent using an ultra-gentle cycle at a maximum of 30°C
- leave to dry away from direct heat
- **NAPOVEDA**
- Handwash with mild detergent at a maximum of 30°C (we recommend using a sponge)
- leave to dry away from direct heat

Dispose of this product and its parts responsibly

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PUTTING ON THE APPLIANCE

- 1 Open the orthosis and temporarily fasten the Velcro® straps on the side (fig. A).
- 2 Put the orthosis in place with the loop on the upper side.
- 3 Insert the central Velcro® strap and fasten well the central Velcro® closure (fig. B). Make sure the closure is centered.
- 4 Tighten the orthosis by fastening the straps on the front Velcro® band (fig. C).
- 5 ONLY on models REF.1091: first tighten the lower straps, then the top ones.

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