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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM**

Elastic orthosis with rigid and thermoformable pads to stabilize the lumbar rachis and relieve pressure

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF.90090	REF.90091*	REF.90092**			
Size	XS	S	M	L	XL	XXL
Pelvic circumf. cm	70/80	80/90	90/100	100/110	110/120	120/130
Front height cm	19	19	19	19	19	19
Back height cm	28	28	28	28	28	28
*Front height cm	23	23	23	23	23	23
*Back height cm	33	33	33	33	33	33
**Front height cm	23	23	23	23	23	23
**Back height cm	55	55	55	55	55	55

Colour: grey

OPTION

REF. 90090 9ST012/L Optional rigid stays cm 22

REF. 90091 9ST012/L Optional rigid stays cm 26

REF. 90092 9ST012/L Optional rigid stays cm 46

MAINTENANCE

 Do not bleach  No chemical cleaning

 Do not iron  Do not tumble-dry

 Washing instructions: hand wash at max. 30° C with neutral soap.

Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS REF. 90090/90091

- Lumbago
- Sciatic lumbalgia and lumbar cruralgia
- Spondylarthrosis and lumbar discopathies
- Slight traumas to the lumbosacral column and paravertebral lumbar contractures
- Post-operative recovery
- REF.90091 only: After-effects of stable compression fractures of lumbar vertebrae
- REF.90091 only: Treatment of fractures to transverse processes

CHARACTERISTICS AND MATERIALS REF. 90090/90091

- Structure in CarbonTrace® fabric, elastic, breathable, ladder-proof, antistatic and bacteriostatic
- Elastic straps for variable and personalised support
- Velcro® front fastening
- 2 front and back thermoformable pads, on medical prescription
- Optional rigid stays, mouldable on medical prescription; to insert in the 4 linear back pockets
- The abdominal section is reinforced by spiral stays fixed in an oblique position
- Individually adaptable to the patient's pathology

INDICATIONS REF. 90092

- Treatment of back pain caused by diffuse spinal osteoarthritis
- Rheumatism
- Post-surgical protection of vertebral fractures treated surgically with instruments and for dorsal herniotomy

CHARACTERISTICS AND MATERIALS REF. 90092

- Padded shoulder straps
- Structure in CarbonTrace fabric, elastic, breathable, ladder proof, antistatic and bacteriostatic
- Elastic straps for variable and personalised support
- Velcro® closure on the front
- The abdominal section is reinforced by spiral stays fixed in an oblique position
- Dorso-lumbar part supported by adjustable stays, integrated into the belt
- 2 front and back thermoformable pads, on medical prescription
- Optional rigid stays, mouldable on medical prescription; to insert in the 4 linear back pockets
- Individually adaptable to the patient's pathology

CONTRAINDICATIONS

Currently no known

PUTTING ON THE APPLIANCE REF. 90090/90091

- 1 Open the corset and temporarily fasten the Velcro® straps on the side (fig. A).
- 2 Put the corset in place with the loop on the upper side.
- 3 Insert your thumb in the loop and fasten well the central Velcro® closure (fig. B). Make sure the closure is centered.
- 4 Tighten the corset by fastening the straps on the front Velcro® band (fig. C). ONLY on models REF.90091: first tighten the lower straps, then the top ones.
- 5 If the use of the lumbar pad (REF.1078/P) is prescribed, fix the pad on the inside to the back Velcro® before putting the corset on.

PUTTING ON THE APPLIANCE REF. 90092

- 1 Open the corset and temporarily fasten the Velcro® straps on the side (fig. D). Put the corset in place with the loop on the upper side.
- 2 Insert your thumb in the loop and fasten well the central Velcro® closure (fig. E).
- 3 Make sure the closure is centered. If necessary reposition the corset (fig. F).
- 4 Tighten the corset by fastening the lower straps on the front Velcro® band (fig. G) first.
- 5 Then tighten the upper straps and fasten them on the front Velcro® band (fig. H).
- 6 Place the upper tensioning straps over the shoulders, then feed under the arms and cross over at the back. Make sure they go through the loops (fig. I). Then secure them to the front Velcro® band (fig. L).

NOTE:

To apply the rear rigid pad (heat-mouldable):

- apply the adhesive padding to the plastic part (concave side-without any writing, fig. M)
- apply the adhesive Velcro® fasteners centrally on the opposite side to the padding (fig. N)
- attach the pad with Velcro® to the internal dorsal part of the orthosis

To apply the front rigid pad (heat-mouldable):

- apply the adhesive padding to the plastic part (concave side-without any writing)
- attach the pad with Velcro® to the front closure of the orthosis

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