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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

## Low / high elastic orthosis with rigid and thermoformable pads to stabilize the lumbar rachis and relieve pressure

### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

### PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields. Do not apply in direct contact with open wounds.

### WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never after the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

### SELECTION/SIZE

Item	REF.33170					
Size	XS	S	M	L	XL	XXL
Circ. waist cm	50/60	60/75	75/90	90/105	105/120	120/135
Front height cm	20	20	20	20	20	20
Back height cm	26	26	26	26	26	26

Colour: grey

OPTION: REF.9ST012/L optional rigid stays cm 20/22


Item	REF.33172					
Size	XS	S	M	L	XL	XXL
Circ. waist cm	50/60	60/75	75/90	90/105	105/120	120/135
Front height cm	22	22	22	22	22	22
Back height cm	32	32	32	32	32	32

Colour: grey

OPTION: REF.9ST012/L optional rigid stays cm 24/26

### MAINTENANCE

-  Do not bleach
-  No chemical cleaning
-  Do not iron
-  Do not tumble-dry

 Washing instructions: hand wash at max. 30° C with neutral soap.

Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

### INDICATIONS

- Lumbago
- Lumbar sciatica and cruralgia
- Spondyloarthrosis and lumbar discopathy
- Slight traumas of the lumbosacral column and paravertebral lumbar contractures
- Postoperative recovery
- LumboSkill72 only: Results of vertebral compressions of the lumbar spine
- LumboSkill72 only: Treatment of fracture of transverse processes

### KONTRAINDICATIONS

Currently no known

### CHARACTERISTICS AND MATERIALS

- The graduated design of the elastic compression fabric on the front and on the sides provides more compression in the middle and a looser fit around the edges; made of polyester and elastomer
- The back is made from BreathSkill fabric:
  - Soft and breathable
  - 3D weaving
  - special laminated covering for extra resistance
- Crossed strengthening straps provide variable and customisable compression
- The abdominal section is reinforced with angled spiral stays
- Front Velcro® fastening system
- Adapts to the pathology of the patient
- 2 front and back thermoformable pads, on medical prescription
- Optional rigid stays, mouldable on medical prescription; to insert in the 4 linear back pockets

### PUTTING ON THE APPLIANCE

- 1 Open the orthosis and temporarily fix the Velcro® straps at the lateral straps
- 2 Open the orthosis and position with the orange label facing upwards.
- 3 Fit the orthosis and fix using the middle Velcro® closure (fig. A).
- 4 Tighten the orthosis by fixing the Velcro® straps on the front Velcro® strap (fig. B). ONLY high model (LumboSkill72, REF.33172): First pull tight the lower and then the upper Velcro® straps (fig. C).
- 5 Ensure that the orthosis is exactly in the centre of the body, both on the front as well as on the back (fig. D).

### NOTE:

To apply the rear rigid pad (heat-mouldable):

- apply the adhesive padding to the plastic part (concave side-without any writing, fig.E)
- apply the adhesive Velcro® fasteners centrally on the opposite side to the padding (fig.F)
- attach the pad with Velcro® to the internal dorsal part of the orthosis

To apply the front rigid pad (heat-mouldable):

- apply the adhesive padding to the plastic part (concave side-without any writing)
- attach the pad with Velcro® to the front closure of the orthosis

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Velcro® è un marchio di fabbrica registrato di Velcro Industries B.V.