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 **DESIGNED AND DEVELOPED IN SWITZERLAND**

 **ORTHOSERVICE RO+TEN**
Take care feel better

 **ORTHOSERVICE RO+TEN**
Take care feel better

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EN – The descriptions and images in this document are for illustrative and commercial purposes only. Orthoservice reserves the right to make modifications without prior notice based on their needs.

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info@orthoservice.com · www.orthoservice.com

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Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de

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 EN
 +
 DE
 +
 IT
 +
 FR
 +
 ES
 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Padded orthosis for stabilising and supporting the dorso-lumbar rachis, with rigid and thermoformable pads

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation.

In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

MAINTENANCE

 Do not bleach
 No chemical cleaning

 Do not iron
 Do not tumble-dry

 Washing instructions: hand wash at max. 30° C with neutral soap.

Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Back pain due to arthrosis or kyphosis
- After-effects of fractures of dorsal vertebrae and/ or upper lumbar vertebrae
- After-effects of stable compression fractures of the dorsal spine due to osteoporosis or metastasis
- Other indications on medical prescription

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Base fabric:
 - 98% reinforced anaelastic cotton for optimal containment
 - Breathable and natural for best comfort
 - Exclusive design
- Anatomical design shaped for the legs. Available in two sizes:
 - Regular fit (REF. 1335)
 - Tailored fit (REF. 1336)
- Traction system:
 - Reinforced elastic straps for targeted pressure and perfect stabilization of the lumbar rachis, latex free
 - Pulley system for easy fastening
 - Ferrules in a soft polymeric material to facilitate grip and removal of the brace
- Front Velcro® fastening system, with comfortable fast for easy closing. Shaped for maximum comfort.
 - Padded underarm straps
 - 2 rigid front and back thermoformable pads, on medical prescription
 - Optional rigid stays, moldable on medical prescription, to be inserted in the 4 linear back pockets.
- Adapts to the pathology of the patient.

PUTTING ON THE APPLIANCE

- Open up the orthosis and provisionally affix the Velcro® straps to the strips on the side (fig. A). Position the support in such a way that the inside label faces upwards.
- Put on the orthosis by sliding the arms through the shoulder straps as if it were a backpack (fig. B).
- Use the middle Velcro® closure to fix in place securely (fig. C). Ensure that the closure is precisely centered.
- Slowly tighten up the support by securing the lower Velcro® straps to the front strips (fig. D).
- Then secure the upper straps to the front Velcro® strips (fig. E).
- Grasp the straps and pull towards the front (fig. F). Attach them to the Velcro® on the strip at the front (fig. G).

To apply the rear rigid pad (heat-mouldable):

- apply the adhesive padding to the plastic part (concave side-without any writing, fig.H)

- apply the adhesive Velcro® to the sides of the central rib (fig.I)

- attach the pad with Velcro® to the internal dorsal part of the orthosis (fig.L)

To apply the front rigid pad (heat-mouldable):

- apply the adhesive padding to the plastic part (concave side-without any writing)

- attach the pad with Velcro® to the front closure of the orthosis

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