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
 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Orthosis with anterior-posterior traction for stabilising and supporting the dorso-lumbar rachis, with a rigid frame

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient.

To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation.

The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation.

In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF. 1370					
Size	XS	S	M	L	XL	XXL
Pelvic circ. cm	75/85	85/95	95/105	105/115	115/125	120/135
Front height cm	21	21	21	21	21	21
Back height cm	55	55	55	55	55	55
Colour	beige					

MAINTENANCE

 Do not bleach  No chemical cleaning
 Do not iron  Do not tumble-dry
 Washing instructions: hand wash at max. 30° C with neutral soap (a sponge is recommended); do not dry in the vicinity of heat sources (or, better, dry with a cloth).

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Back pain caused by arthrosis or kyphosis
- Recovery from fractures of the thoracic and/or upper lumbar vertebrae
- Management of thoracic vertebral compression fractures caused by osteoporosis or metastasis

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Dorsal traction: an innovative reverse system (anterior-posterior) for optimal support
- Lumbar traction: stretch straps provide targeted pressure and perfect stabilisation
- Front Velcro® fastening
- Shoulder straps with underarm padding
- Adapts to the pathology of the patient
- Cotton fabric with minimal elastane content
- Bonded two-layer MTP shoulder straps
- Rigid reinforced stays provided, mouldable on medical prescription; to insert in the 4 linear back pockets

PUTTING ON THE APPLIANCE

- 1 Undo the corset, attach the Velcro® straps to the side strips for now (fig. A), and attach the straps on the shoulder straps to the white Velcro® markers (fig. B) for the time being.
- 2 Put on the corset by sliding the arms through the shoulder straps as if it were a backpack (fig. C).
- 3 Close the corset at the front by attaching the central Velcro® fastening securely (fig. D). Check that the front fastening is properly centred.
- 4 Start to tighten the corset by attaching the lower straps to the Velcro® strip at the front (fig. E).
- 5 Next attach the top straps to the Velcro® strip at the front (fig. F).
- 6 Take hold of the straps on the shoulder straps (which had previously been attached to the white Velcro® markers) and pull them towards the front (fig.G). Attach them to the strip at the front using the Velcro® provided (fig. H).
- 7 **[ONLY FOR THE FIRST FITTING]** If necessary, adjust the length of the shoulder straps: remove the Velcro® tips (fig. I), cut the straps to the correct length (fig. L) and reattach the Velcro® tips.

To ensure optimal stabilisation, mould the rigid stays provided and slide them into the pockets provided in the back.

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