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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Low/high orthosis with easy-pull tensioning straps for stabilising and supporting the lumbar spine,with rigid thermoformable pad

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF.90070				
Size	S	M	L	XL	XXL
Pelvic circumf. cm	70/85	85/100	100/115	115/130	130/145
Front height cm	16	16	16	16	16
Back height cm	26	26	26	26	26
Colore: grey					

Item	REF.90072				
Size	S	M	L	XL	XXL
Pelvic circumf. cm	70/85	85/100	100/115	115/130	130/145
Front height cm	20	20	20	20	20
Back height cm	31	31	31	31	31
Colore: grey					

MAINTENANCE

 Do not bleach
 No chemical cleaning

 Do not iron
 Do not tumble-dry

 Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources. Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Lumbago
- Sciatic lumbalgia and lumbar cruralgia
- Spondylarthrosis and lumbar discopathies
- Slight traumas to the lumbosacral column and paravertebral lumbar contractures
- Post-operative recovery
- After-effects of stable compression fractures of lumbar vertebrae
- Treatment of fractures to transverse processes
- Lumbopull72 only: Degenerative lumbar scoliosis
- Lumbopull72 only: Spondylosis with or without vertebral listhesis

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Compound fabric (3D fabric on the inside, special fabric on the outside): high rigidity to ensure the greatest possible stabilisation; breathable for the greatest possible comfort when wearing
- System of serially arranged pulleys and tensioning straps (exclusive) in order to achieve optimum stabilisation and tensioning with minimum application of strength. The closure multiplies the force exerted by the patient (x 4) and ensures the even distribution of force along with optimum adaptability to the individual anatomy
- Elastic rear insert (not structural) in order to prevent the back coming into contact with the straps and to apply various optional items using Velcro® fastener
- Rigid paravertebral reinforcements made of polymer material to ensure stability and rigidity, thermoformable
- Anatomically shaped and easy to put on
- Large adhesive surface for maximum setting options
- Terminals for strap traction with anatomical and facilitated grip
- Low rigid pad, thermoformable, upon medical prescription
- Optional:
 - high rigid pad (Pad20)
 - back support (DorsoPull74).

PUTTING ON THE APPLIANCE FOR THE PATIENT

- (First time only) If prescribed by the doctor, attach the rigid lumbar pad to the inside of the middle rear of the corset using the Velcro® closure (fig.A1); it is possible to thermoform it to adapt it to the specific needs of the patient (fig.A2).
- Open the corset (make sure that the tensioning straps are pulled out to their full length) and secure the Velcro® end pieces temporarily to the sides (fi g. B).
- Position the corset in such a way that the label on the inside points upwards and, when putting on, use the centre Velcro® closure to fasten securely (fig. C).

Note: Pull the corset tight.
- Ensure that the central Velcro® closure is positioned in the middle.
- Tension the corset by pulling on the strap-pulley system using the designated Velcro® closure ends (fig. D) and secure these to the front adhesive surface using the Velcro® (fig. E).

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Velcro® è un marchio di fabbrica registrato di Velcro Industries B.V.