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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Padded shoulder back splint for the orthopedic Lumbopull orthosis

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device. Do not apply in direct contact with open wounds

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF.90074				
Size	XS	S	M	L	
Distance T1-L5 cm	35/41	41/47	47/54	54/61	
Colour	grey				

MAINTENANCE

- Do not bleach
- No chemical cleaning
- Do not iron
- Do not tumble-dry
- Washing instructions: hand wash at max. 30° C with neutral soap.

Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Back pain due to arthrosis or kyphosis
- After-effects of fractures of dorsal vertebrae and/or upper lumbar vertebrae
- After-Effects of stable compression fractures of the dorsal spine due to osteoporosis or metastasis

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- When combined to Lumbopull 70 o 72, it becomes a padded shoulder orthopedic orthosis corset with shoulder straps
- Self-shaping dorsal stays included
- There is also the option of inserting a formable, stable internal frame which is made of aluminium alloy and radiolucent
- New length adjusting shoulder straps with an elastic panel at the back enable a secure fit and the brace to be worn like a backpack
- Two-fold compound padding made of PE foam

ADAPTING FOR THE DOCTOR/TECHNICIAN

- 1 Remove the rigid inner frame from the rear section (fig. A) adjust it to the patient's anatomy and put it back into the rear section: for the adjustment, use, for example, the edge of a table or other surface for leverage if necessary (or use an upset) (fig. B).
- 2 Attach the rear brace using the Velcro® closure to the inside of the centre section of the corset (fig. C). Position the point at which the straps cross at the desired height using the Velcro (fig. D1) and pass the free ends of the dorsal-lumbar straps through the loops provided on the corset (fig. D2).
- 3 Open the corset (make sure that the tensioning straps are pulled out to their full length) and secure the Velcro® end pieces of the corset and back brace temporarily to the sides (fig. E).
- 4 The shoulder straps of the brace are put on in a similar way to the carry straps of a rucksack (fig. F).
- 5 Put on the corset and use the middle Velcro® closure to fasten securely (fig. G).
Note: Pull the corset tight.
- 6 Ensure that the central Velcro® closure is positioned in the middle.
- 7 Tension the corset by pulling on the strap-pulley system using the designated Velcro® closure ends and secure these to the front adhesive surface using Velcro® (fig. H).
- 8 Adjust the length of the shoulder straps (at the top) using the Velcro® fastenings provided (fig. I). If necessary, adjust the position of the padded shoulder cuff in such a way that the tube-shaped section is positioned in the armpit. These padded cuffs have a small adhesive surface on the upper section that enables them to be secured into the optimum position.
- 9 Position the end pieces of the dorsolumbar tensioning straps on the stomach pad, check the length of the straps for tension and cut off any protruding parts (fig. L).

PUTTING ON THE APPLIANCE FOR THE PATIENT

- 1 Open the corset (make sure that the tensioning straps are pulled out to their full length) and secure the Velcro® end pieces of the corset and back brace temporarily to the sides (fig. M).
- 2 The shoulder straps of the brace are put on in a similar way to the carry straps of a rucksack (fig. N).
- 3 Put on the corset and use the middle Velcro® closure to fasten securely (fig. O).
Note: Pull the corset tight.
- 4 Ensure that the central Velcro® closure is positioned in the middle.
- 5 Tension the corset by pulling on the strap-pulley system using the designated Velcro® closure ends and secure these to the front adhesive surface using Velcro® (fig. P).
- 6 Close the shoulder straps. To do this, pull the straps outwards and secure to the adhesive surface on the front (fig. Q).

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