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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Semi-rigid corset with lumbar pelotte or Chairback-frame

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Size	REF. T1030R-T1031R / T1030E-T1031E / T1030C-T1031C					
Taglie	S	M	L	XL	XXL	XXXL
Belt circ. cm	79/90	90/103	103/115	115/128	128/144	144/160
Height (REF. T1030)	13/21	13/21	15/22	18/23	20/24	23/25
Height (REF. T1031)	19/25	20/25	20/25	20/25	23/26	23/26
Colour	black					

Option: abdominal pelotte, easily fixable using the Velcro® strap

MAINTENANCE

- Do not bleach
- Do not iron
- No chemical cleaning
- Do not tumble-dry
- Washing instructions: hand wash at max. 30°C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

SCOPE OF DELIVERY

- THE LUMBOSAFE INCLUDES: 1 x belt (low version, REF. T1030 or high version REF. T1031)
- 1-2 x small abdominal/lumbar pelottes, „Regular“ (additional pelotte only included in REF. T1030R • T1031R)
- 1 x large lumbar pelotte „Extra“ (only for REF. T1030C • T1031E)
- 1 x „Chairback“ pelotte (only for REF. T1030C • T1031C)

INDICATIONS

- Lumbalgia
- Sciatic lumbalgia and lumbar cruralgia
- Spondylarthrosis and lumbar discopathies
- Slight traumas to the lumbosacral column and paravertebral lumbar contractures
- Post-operation period
- Results of vertebral lumbar relaxations
- Treatment of fracture of transverse processes

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Innovative tensioning system with pulley technology for optimum stabilisation and for simple application with a minimum of force required
- Firm, form-stable and breathable fabric for safe hold, with border of skin-friendly material. Particularly light
- The large Velcro® closure surface enables the tensioning strap to be simply secured in variable pressure settings
- Reinforced paravertebral area
- Stable, easily secured, padded back pelotte made from hard plastic which can be quickly replaced with a just a few steps according to therapy form required
- Belt system of two different heights

ADAPTING FOR THE DOCTOR / TECHNICIAN

- If required, attached the prescribed accessories:
 - only for REF. T1030R • T1031R: Attach the half-height lumbar pelotte on the rear part of the belt by feeding the elasticated tapes through the cut-outs designed for this purpose in the plastic component of the pelotte (fig. A)
 - only for REF. T1030E • T1031E: Attach the tall lumbar pelotte on the rear part of the belt by feeding the elasticated tapes through the cut-outs designed for this purpose in the plastic component of the pelotte (fig. B)
 - Only for REF. T1030C • T1031C: Feed the belt into the two Velcro® loops of the chairback pelotte intended for this purpose (fig. C).
- Open the belt and provisionally attach the Velcro® closure element to the side of the belt.
- Position the belt in such a way that the labels „right/left“ are on the correct sides and are pointing upwards (fig. D).
- Put on the belt and secure using the central Velcro® closure (fig. E).
- Optimise the tension of the corsets using the end piece (fig. F) and secure via the Velcro® closure on the appropriate Velcro® surface (fig. G).
- Make sure that the reinforced section at the back are positioned centrally / symmetrically to the spine.

PUTTING ON THE APPLIANCE FOR THE PATIENT

- Open the belt and provisionally attach the Velcro® closure element to the side of the belt.
- Position the belt in such a way that the labels „right/left“ are on the correct sides and are pointing upwards (fig. H).
- Put on the belt and secure using the central Velcro® closure (fig. I).
- Optimise the tension of the corsets using the end piece (fig. L) and secure via the Velcro® closure on the appropriate Velcro® surface (fig. M).
- Make sure that the reinforced section at the back are positioned centrally / symmetrically to the spine.

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