Tutore dorso-lombare con busto SAT

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SINGLE PATIENT
MULTIPLE USE

 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Dorso-lumbar brace with SAT orthosis

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.2709					
Size	XS	S	M	L	XL	XXL
Pelvis circumference cm	70/80	80/90	90/100	100/110	110/120	120/130
Height at back cm	40	40	45	45	50	50
Colour	beige					

MAINTENANCE

 Do not bleach  No chemical cleaning
 Do not iron  Do not tumble-dry
 Washing instructions: hand wash at max. 30° C with neutral soap (a sponge is recommended); do not dry in the vicinity of heat sources.
Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Recovery from fractures of the dorso-lumbar spine
- Osteomalacia
- Dorso-lumbar vertebrae compression fractures caused by osteoporosis
- Vertebral compression fractures caused by metastases
- Spondylitis
- Spondylodiscitis
- Roundback and hyperkyphosis

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- A semi-rigid breathable fabric orthosis (3% elastane).
- Aluminium alloy back frame (removable).
- The traction system with wide elasticated straps provides targeted pressure and optimal stabilisation of the lumbar spine.
- The shoulder straps are adjustable and have moveable padding.
- A rigid pelvic band with non-stretch straps.
- Front Velcro® fastening.

PUTTING ON THE APPLIANCE ADAPTING FOR THE DOCTOR/TECHNICIAN

- 1 If necessary, open the lower Velcro® pocket at the back and remove the metal frame from the fabric orthosis (fig. A). Mould the frame to the desired shape and slide it back into the orthosis.
- 2 Undo the orthosis and attach the Velcro® straps of the elastic straps to the side strips for the time being (fig. B). Attach the ends of the shoulder straps to the Velcro® areas provided on the pelvic band (fig. C).
- 3 Slide the arms through the shoulders straps and put on the orthosis (fig. D).
- 4 Close the orthosis by inserting a thumb in the loop and securely fastening the central Velcro® fastener (fig. E). Check that the front fastening is properly centred.
- 5 Start to tighten the orthosis by fastening the Velcro® on the elastic straps to the Velcro® strip at the front (fig. F).
- 6 Next pass the belt of the pelvic band through the ring provided, pull it tight and fasten it back on itself (fig. G).
- 7 Adjust the length of shoulder straps using the loops provided (fig. H). Position the shoulder straps so they protect the shoulders and armpits (fig. I).
- 8 Fasten the ends of the shoulder straps by pulling them outwards (fig. L) and then fastening them with Velcro® to the front.

PUTTING ON THE APPLIANCE FOR THE PATIENT

- 1 Undo the orthosis and attach the Velcro® straps of the elastic straps to the side strips for the time being (fig. M). Attach the ends of the shoulder straps to the Velcro® areas provided on the pelvic band (fig. N).
- 2 Slide the arms through the shoulders straps and put on the orthosis (fig. O).
- 3 Close the orthosis by inserting a thumb in the loop and securely fastening the central Velcro® fastener (fig. P). Check that the front fastening is properly centred.
- 4 Start to tighten the orthosis by fastening the Velcro® on the elastic straps to the Velcro® strip at the front (fig. Q).
- 5 Next pass the belt of the pelvic band through the ring provided, pull it tight and fasten it back on itself (fig. R).
- 6 Fasten the ends of the shoulder straps by pulling them outwards (fig. S) and then fastening them with Velcro® to the front.

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