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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Spinal brace for osteoporosis

KNOWLEDGE OF THE USER: The ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

IMPORTANT INFORMATION: Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Use the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin can cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF.33049				
Size	XS	S	M	L	
Distance T1-L5 cm	35/41	41/47	47/54	54/61	
Pelvic circ. cm	70/125	70/125	70/125	70/125	
Colour	grey				

REF. 33049/P optional inflatable pad

MAINTENANCE

- Do not bleach
- No chemical cleaning
- Do not iron
- Do not tumble-dry
- Washing instructions: Pull out the metal frame, see points 1-a, b, c. Wash by hand in lukewarm water with neutral soap, using thoroughy. Leave to dry away from heat sources. Reinsertion of the metal frame: See points 1-f, g, h. Do not dispose of the product or any of its components into the environment.

- Dorsolumbalgia after healing of stable compression fractures of dorsolumbar vertebrae
- Stable compression fractures of dorsolumbar vertebrae caused by osteoporosis or metastatic tumors
- Hyperkyphosis of the thoracic spine, post-traumatic or caused by arthrosis
- Juvenile osteochondrosis

CHARACTERISTICS AND MATERIALS

- Malleable radiotranslucent aluminum alloy inner frame
- Patented strap system
 - Fastens around the abdomen with a corset
 - A unique tensioning system – all the rear straps and the shoulder straps can be tightened by pulling the pads
 - A special system magnifies the force exerted to double the tension of the shoulder straps – maximum effectiveness with minimum effort
- Bonded straps with low-friction strap
- Dorsal-abdominal straps with attached-rixtion to the metal frame with self aligning buckles and press studs
- The abdominal corset ensures greater comfort and stability:
 - Velcro® fastening
 - Adjustable length
 - elastic stretch MTP section and technical stretch fabric
 - Biodynamic, stretchy dorsal insert
 - Padded shoulder straps with the option of pre-setting their positions
 - PE padding bonded with velour and technical fabric with bonded edges

ADAPTING FOR THE DOCTOR/TECHNICIAN

- Remove the metal frame from the padding:
 - Unfasten the four press studs (fig. A) that attach the pad to the metal frame.
 - Unfasten the top buckles (fig. B) and push them so they are released from the pin (N.B.: the buckles have a pin to prevent them from coming out of the fabric).
 - Slide out the metal part (frame) from the padding (fig. C).
 - Rest the metal frame on the patient's back and check if any points need to be shaped (fig. D).
 - Fit the brace in accordance with the medical instructions: using a surface as leverage can make it easier to shape the brace, e.g. the edge of a table (fig. E).
 - Slide the metal part back into place in the padding (fig. F) and ensure that the pins are at the top.
 - Fasten the two lower press studs (fig. G1).
 - Fasten the two central press studs (fig. G2) – N.B.: there are two possible positions for the central attachment to suit the size of the patient's shoulders; use the indentations on the outer padding as a guide when aligning the press studs.
 - Fasten the top buckles (fig. H), ensuring that they fit onto the pins, and pull them until you hear a click.
 - The brace is now ready to be fitted on the patient; please note that the shoulder straps and the straps are fully extended at this point.
- To shorten the corset:
 - Undo the Y-shaped fastenings and remove the abdominal corset from the dorsal part (fig. K1).
 - Measure the circumference of the patient's pelvis.
 - Cut off an excess from the velour part using the following table as a guide: **TAKE CARE not to cut too much.**

Pelvis circumference. cm	70	80	90	100	110
Minimum amount to cut off each of the two sides cm	25	20	15	10	5

- Reattach the corset to the dorsal part using the Velcro® fastenings (fig. K1).
- To put on the brace, slide the patient's arms through the shoulder straps (fig. J).
- Fasten the abdominal pads with Velcro® so that they are centred on the abdomen in a slightly staggered position (fig. I):
 - Check that the rear part of the brace is centred over the spine.
 - Check that the brace is further shortened as described in point 2.
- Adjust the shoulder straps using the tip with the Y-tab. N.B.: it is very important that the ring which attaches to the dorsal-lumbar strap is centred on the median line of the body (fig. 1). If a strap is too long, cut to the correct length (fig. M).
- Use the Velcro® fastenings provided to adjust the positioning of the pads on the shoulder straps so that the tubular part remains under the armpit (fig. N).
- Position the pads of the dorsal-lumbar straps on the abdominal pads (fig. O), pass the straps inside these until the desired degree of tension is achieved (fig. P) and slide out the corresponding card to attach the straps (fig. Q). Check that the tension is correct by positioning the pads in the centre of the abdominal pad and adjust further if required. Only cut off any excess when this procedure has been completed (fig. R).
- Check that the straps pass through the corresponding holes on the corset (fig. S) so that it is easier to don in the future.
- As prescribed, attach the OPTIONAL inflatable pad (REF. 33049PA) to the internal Velcro® fastening at the back (fig. T) and ensure it is correctly positioned and inflated as required to increase dorsal pressure. N.B. – do not exert too much pressure on the spine.

PUTTING ON THE APPLIANCE FOR THE PATIENT

- Fitting the orthosis**
 - Fit the orthosis by inserting the arms in the shoulder straps (fig. 1).
 - Loose the abdominal pelottes using the Velcro® closure. Ensure that they are centred (fig. 2).
 - Loosen the pelottes of the tensioning straps of the back and lumbar area (pre-set at the rear ends of the abdominal tensioning straps), pull outwards (fig. 3) and fix on the front closure using the Velcro® closure (fig. 4).
- Taking off the orthosis**
 - Loosen the small pelottes of the tensioning straps of the back and lumbar areas (fig.5) and fold them provisionally on themselves.
 - Loosen the shoulder straps pulling them forward (fig. 6).
 - Loosen the abdominal pelotte (fig. 7).
 - Take off the orthosis (fig. 8).

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MARCHI COMMERCIALI DEI MATERIALI: Velcro® è un marchio di fabbrica registrato di Velcro Industries B.V.