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Table with columns for Art.-Nummer, Größe, Abstand T1-L5 cm, and Farbe. Includes size chart for XS, S, M, L and a small image of the brace.

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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM
Spinal brace for osteoporosis

DECLARATION OF CONFORMITY
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.
TRADEMARKS OF MATERIALS
Velcro® is a registered trademark of Velcro Industries B.V.
SAFETY PRECAUTIONS
We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.
WARNINGS
It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.
SELECTION/SIZE

Item REF.T1049 XS S M L
Size 35/41 41/47 47/54 54/61
Distance T1-L5 cm
Colour blue

Abdominal and dorsal tensioning straps, made from comfortable belt material
MAINTENANCE
Do not bleach
Do not iron
Do not tumble-dry
No chemical cleaning
Do not machine wash
Washing instructions: Pull out the metal frame, see points 1-a, b, c. Wash by hand in lukewarm water with neutral soap, rinse thoroughly. Leave to dry away from heat sources. Reinsertion of the metal frame: See points 1-f, g, h. Do not dispose of the product or any of its components into the environment.
SCOPE OF DELIVERY
A Shoulder straps
B Upper buckle with automatic blocking anchorage
C Shoulder strap padding
D Padded body with embedded metal frame
E Tensioning straps of the back and lumbar area
F Snap fasteners
G Abdominal tensioning strap
H Small pelotte of the tensioning straps of the back and lumbar area
I Abdominal pelotte
INDICATIONS
Dorsolumbalgia after healing of stable compression fractures of dorsolumbar vertebrae
Stable compression fractures of dorsolumbar vertebrae caused by osteoporosis or metastatic tumours
Hyperkyphosis of the thoracic spine, posttraumatic or caused by arthrosis
Juvenile osteochondrosis
CHARACTERISTICS AND MATERIALS
Padded shoulder straps with position settings which can be set prior to use
Velcro® closure
Elastic insert with biodynamic effect
Mouldable inner frame in aluminium alloy, X-ray permeable
Abdominal and dorsal tensioning straps, made from comfortable belt material
Buckles with self-aligning mounting points
State-of-the-art system of tensioning straps, protected by patent
Two-layer PE padding with suede and mesh fabric cover
ADAPTING FOR THE DOCTOR/TECHNICIAN
1 Remove the metal frame from the padding:
a. Open the four snap fasteners (fig. A), which fix the padding in the metal frame.
b. Unhook the upper buckle (fig. B), by pressing on the buckles so that they can be removed from the pin (note: The buckles are equipped with a stop mechanism so that they do not slip out of the fabric).
c. Pull the metal frame out of the padding (fig. C).
d. Place the metal frame at the back of the patient to determine, what modifications have to be carried out (fig. D).
e. Proceed then with the modelling in accordance with the medical indications: Use e.g. the edge of a table or a surface that can be used to generate a leverage effect as well as bending irons (fig. E).
f. Reinsert the metal frame into the padding (fig. F) and ensure that the pins are facing upwards.
g. Re-close the four snap fasteners (fig. G).
h. Hook the upper buckles in the pin (fig. H) and tighten until you hear a click. The orthosis is now ready for fitting on the patient.
2 Tighten the left and right abdominal tensioning straps, adjust the length and cut any protruding parts (ideally the straps should not jut out from the pelottes after the final fixing), pull at the Y-shaped ends until the correct tension is achieved and fix them (fig. L).
3 Set the shoulder straps using the Y-shaped ends. Attention: It is very important that the connection rings of the tensioning straps of the back and lumbar area are centred on the median line of the body (fig. M). Cut any protruding parts (fig. N).
4 Position the pads in such a way that the tube-shaped part is located in the armpit. The pads have a small Velcro® closure in the upper part which can be used to fix the pads in the ideal position (fig. O).
5 Position the small pelottes of the tensioning straps of the back and lumbar area on the abdominal pelottes (fig. P), check the length of the tensioning straps when tensioned and cut the protruding part (fig. R). Check the optimum tension by positioning the small pelottes in the centre of the abdominal pelotte and, if required, carry out some fine adjustment (fig. S).
6 Check the overall fit of the orthosis and tensioning straps of the shoulder strap-tensioning strap system by checking the strap lengths (fig. T). Check that shoulder straps and tensioning straps have the right length, and cut, if required.

Table with columns for Item, Size, Distance T1-L5 cm, and Colour. Includes size chart for XS, S, M, L and a small image of the brace.

Abdominal and dorsal tensioning straps, made from comfortable belt material
MAINTENANCE
Do not bleach
Do not iron
Do not tumble-dry
No chemical cleaning
Do not machine wash
Washing instructions: Pull out the metal frame, see points 1-a, b, c. Wash by hand in lukewarm water with neutral soap, rinse thoroughly. Leave to dry away from heat sources. Reinsertion of the metal frame: See points 1-f, g, h. Do not dispose of the product or any of its components into the environment.
SCOPE OF DELIVERY
A Shoulder straps
B Upper buckle with automatic blocking anchorage
C Shoulder strap padding
D Padded body with embedded metal frame
E Tensioning straps of the back and lumbar area
F Snap fasteners
G Abdominal tensioning strap
H Small pelotte of the tensioning straps of the back and lumbar area
I Abdominal pelotte
INDICATIONS
Dorsolumbalgia after healing of stable compression fractures of dorsolumbar vertebrae
Stable compression fractures of dorsolumbar vertebrae caused by osteoporosis or metastatic tumours
Hyperkyphosis of the thoracic spine, posttraumatic or caused by arthrosis
Juvenile osteochondrosis
CHARACTERISTICS AND MATERIALS
Padded shoulder straps with position settings which can be set prior to use
Velcro® closure
Elastic insert with biodynamic effect
Mouldable inner frame in aluminium alloy, X-ray permeable
Abdominal and dorsal tensioning straps, made from comfortable belt material
Buckles with self-aligning mounting points
State-of-the-art system of tensioning straps, protected by patent
Two-layer PE padding with suede and mesh fabric cover
ADAPTING FOR THE DOCTOR/TECHNICIAN
1 Remove the metal frame from the padding:
a. Open the four snap fasteners (fig. A), which fix the padding in the metal frame.
b. Unhook the upper buckle (fig. B), by pressing on the buckles so that they can be removed from the pin (note: The buckles are equipped with a stop mechanism so that they do not slip out of the fabric).
c. Pull the metal frame out of the padding (fig. C).
d. Place the metal frame at the back of the patient to determine, what modifications have to be carried out (fig. D).
e. Proceed then with the modelling in accordance with the medical indications: Use e.g. the edge of a table or a surface that can be used to generate a leverage effect as well as bending irons (fig. E).
f. Reinsert the metal frame into the padding (fig. F) and ensure that the pins are facing upwards.
g. Re-close the four snap fasteners (fig. G).
h. Hook the upper buckles in the pin (fig. H) and tighten until you hear a click. The orthosis is now ready for fitting on the patient.
2 Tighten the left and right abdominal tensioning straps, adjust the length and cut any protruding parts (ideally the straps should not jut out from the pelottes after the final fixing), pull at the Y-shaped ends until the correct tension is achieved and fix them (fig. L).
3 Set the shoulder straps using the Y-shaped ends. Attention: It is very important that the connection rings of the tensioning straps of the back and lumbar area are centred on the median line of the body (fig. M). Cut any protruding parts (fig. N).
4 Position the pads in such a way that the tube-shaped part is located in the armpit. The pads have a small Velcro® closure in the upper part which can be used to fix the pads in the ideal position (fig. O).
5 Position the small pelottes of the tensioning straps of the back and lumbar area on the abdominal pelottes (fig. P), check the length of the tensioning straps when tensioned and cut the protruding part (fig. R). Check the optimum tension by positioning the small pelottes in the centre of the abdominal pelotte and, if required, carry out some fine adjustment (fig. S).
6 Check the overall fit of the orthosis and tensioning straps of the shoulder strap-tensioning strap system by checking the strap lengths (fig. T). Check that shoulder straps and tensioning straps have the right length, and cut, if required.

PUTTING ON THE APPLIANCE FOR THE PATIENT
Fitting the orthosis
1 Fit the orthosis by inserting the arms in the shoulder straps (fig. 1).
2 Close the abdominal pelottes using the Velcro® closure. Ensure that they are centred (fig. 2).
3 Loosen the pelottes of the tensioning straps of the back and lumbar area (pre-set at the rear ends of the abdominal tensioning straps), pull outwards (fig 3) and fix on the front closure using the Velcro® closure (fig. 4).
Taking off the orthosis
4 Loosen the small pelottes of the tensioning straps of the back and lumbar areas and fix provisionally at the rear ends of the abdominal tensioning straps (fig. 5).
5 Loosen the abdominal pelotte (fig. 6).
6 Loosen the shoulder straps (fig. 7).
7 Take off the orthosis (fig. 8).

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Table with columns for Code, Mesure, Distance T1-L5 cm, and Couleur. Includes size chart for XS, S, M, L and a small image of the brace.

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