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Currently no known

CHARACTERISTICS AND MATERIALS

- A lightweight aluminium alloy frame without any sharp metal edges for maximum comfort.
- The brace can be adjusted at the joints with screws and at the top with a telescopic mechanism (the screws do not need to be removed).
- The angle of the underarm girdles can be adjusted using the notched system that corresponds with the angle of the underarm pads (7,5° increments).
- An innovative polymer spring enables the pelvic band to be tilted (it can also be secured in position at 0° and at 18° using grub screws).
- Moveable sternal plate.
- A new strap with three-step adjustable length and a clip to hold any excess strap.
- The padded parts are water-repellent and can be immersed in water: PE foam covered with non-toxic synthetic leather.
- The back plate can be worn horizontally and vertically.

PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

3-Point hyperextension brace

DECLARATION OF CONFORMITY
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS
Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS
We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. In doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition, the external label carefully. We recommend not to wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS
It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or of any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic device can only be guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF P1000
Size	XS S MS M ML LS L XLS XL
Pelvic circ. cm	50/60 60/75 75/90 75/90 75/90 90/105 90/105 105/115 105/115
Height cm*	30/37 34/42 34/41 36/46 40/52 39/49 42/54
Colour	grey

* measured from the proximal side of the corpus sterna to the upper edge of the pubic bone (tuberculum pubicum)

MAINTENANCE

- ⚠ Do not bleach ☒ No chemical cleaning
- 🌀 Do not iron ☒ Do not tumble dry
- 🧼 Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Traumatic fractures (T10-L1/L2)
- Osteoporotic or metastatic-related compression fractures of the thoracic spine and the dorsolumbar transition (T10-L1/L2)
- Osteomalacia with fractures
- Dorsolumbar osteoarthritis
- Substitute for plaster corset

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- A lightweight aluminium alloy frame without any sharp metal edges for maximum comfort.
- The brace can be adjusted at the joints with screws and at the top with a telescopic mechanism (the screws do not need to be removed).
- The angle of the underarm girdles can be adjusted using the notched system that corresponds with the angle of the underarm pads (7.5° increments).
- An innovative polymer spring enables the pelvic band to be tilted (it can also be secured in position at 0° and at 18° using grub screws).
- Moveable sternal plate.
- A new strap with three-step adjustable length and a clip to hold any excess strap.
- The padded parts are water-repellent and can be immersed in water: PE foam covered with non-toxic synthetic leather.
- The back plate can be worn horizontally and vertically.

ADAPTING FOR THE DOCTOR/ORTHOPEDIC TECHNICIAN

- 1 Place the hyperextension brace on the front of the patient's chest (preferably with the patient lying down) to check which adjustments are required (fig. A).
- 2 Loosen the screws on the vertical side stays (fig. B1) and on the horizontal bar (fig. B2).
- 3 Adjust the horizontal bar as follows:
 - Adjust the brace to the width of the patient's chest by moving the two overlapping segments (fig. C1);
 - Tighten the screws, checking that the sternal pad is symmetrical with the sagittal plane (fig. C2);
- 4 Next adjust the vertical stays as follows:
 - Lengthen/shorten the stay so that the pelvic band is level with the pubis and the sternal plate is level with the sternum (fig. D1); remove the screws if necessary, and then reposition them.
 - Tighten the screws (fig. D2) and check that the patient can sit down without any discomfort; adjust the vertical stays again if necessary.
 - Check that the two stays are the same length.
- 5 Adjust the angle of the underarm joints (between the vertical stays and the horizontal bar):
 - Remove the screw and the metal cover (fig. E1).
 - Lift the underarm part, turn it to the desired angle, and then reposition it on the notched disc (fig. E2).
 - Replace the metal cover and screws; tighten the screws.
 - Repeat the same procedure on the other side, ensuring that it is adjusted to the same angle.
- 6 Adapt the brace to fit the patient's body exactly:
 - Shape the pelvic band to fit (fig. F1).
 - Shape the underarm pads (fig. F2), ensuring that they fit snugly against the chest area.
 - Shape the segments at the height of the pelvis (fig. F3), adjusting the angle of inclination.
- 7 If specified by the doctor, place the lumbar plate against the spine in the position prescribed:
 - Remove the mushroom head screw and position it at the required height (in one of the brass threaded bushings provided) and tighten (fig. G1).
 - Also move the lever to the desired height using the screws provided (fig. G2).
- 8 If specified by the doctor, the lumbar plate can be positioned vertically (fig. H) by removing the strap from the horizontal keepers and passing it through the vertical keepers.
- 9 Adjust the strap as follows:
 - Open the lever fastening, pass the strap behind the patient's back and insert the hook into the mushroom head screw (fig. I1).
 - Adjust the length of the strap so that it provides the right degree of pressure in the lumbar area: slide the strap through the three-step buckle (fig. I2) and fasten by pulling the end (fig. I3); optimum pressure is achieved when the hyperextension brace stays in the correct position even when the lever fastening is open.
 - If the strap is too long, attach it to size or fasten it open.
 - Click the fastening lever into position (fig. I5).
 - Check that the lumbar plate is centred securely on the spine; fasten it securely to the strap using the keepers provided (fig. I6).

SECURING THE PELVIC BAND (for the doctor/orthopedic technician)

The brace is supplied with a pelvic band that tilts and moves. However, it can be secured using the joints at the side if necessary:

- 1 - Move the pelvic band so that one of the two hinges in the mechanism matches up with the corresponding threaded hole (fig. L) and lock it at 0° or 18°.
- 2 - Insert the grub screw provided and tighten (fig. M).
- 3 - Repeat the same procedure on the opposite side.
- 4 - Check that both mechanisms are secured in the same way.

PUTTING ON THE APPLIANCE FOR THE PATIENT

- 1 Open the lever fastening and unhook the strap from the mushroom head screw.
- 2 Place the hyperextension brace on the chest so the chest side metal plate rests on the centre of the chest and the pelvic band rests on the pubis (fig. N).
- 3 Pass the strap behind the back (fig. O) and hook it to the mushroom head screw (fig. P).
- 4 Secure the hyperextension brace using the fastening lever (fig. Q).

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