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**EN PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

**Cruciform spinal hyperextension orthosis**

**DECLARATION OF CONFORMITY**  
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.  
**TRADEMARKS OF MATERIALS**  
Velcro® is a registered trademark of Velcro Industries B.V.

**SAFETY PRECAUTIONS**  
We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields. Do not apply in direct contact with open wounds. Do not release the tension of the belt abruptly, in order to prevent sudden collapse of the trunk.

**WARNINGS**  
It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

**SELECTION/SIZE**

Code	REF.D100		
Size	S	R	L
Distance sterno pube cm	27/34	34/41	40/53
Waist max circ. cm	90	120	120
Colour	white		

**MAINTENANCE**  
⚠ Do not bleach ⚠ No chemical cleaning  
⚠ Do not iron ⚠ Do not tumble-dry  
⚠ Washing instructions: hand wash at max. 30° C with neutral soap (a sponge is recommended); do not dry in the vicinity of heat sources (or, better, dry with a cloth)  
• Check screws for tight fit regularly  
• Do not dispose of the product or any of its components into the environment.  
• Contact an orthopaedic technician for the replacement of used parts

**INDICATIONS**  
• Extension and immobilisation of the dorso-lumbar spine  
• Post-operative immobilisation  
• Hyperkyphosis and roundback  
• Scheuermann disease  
• Collapsed vertebra caused by osteoporosis  
• Collapsed vertebra caused by metastasis

**CONTRAINDICATIONS**  
Currently no known

**CHARACTERISTICS AND MATERIALS**  
• Lightweight T6 aluminum alloy frame (allows a weight of less than 700 g)  
• Adjustable and movable sternal and pubic pads  
• Lumbar pelotte in 3 sizes: small (horizontal) - regular - long  
• Replaceable EVA padding- waterproof, non toxic and latex-free  
• New PLUS fastening made from light, durable compound material. Simple to open and close, guarantees constant lumbar pressure. With safety fastening  
• Simple telescopic adjustment of frame - no need to remove screws

**ADAPTING FOR THE DOCTOR / TECHNICIAN**

1 Undo the screws in order to regulate the vertical and horizontal position of the hyperextension orthosis (fig. A and B).

2 Have the patient lie on their back and position the hyperextension orthosis on the front of the torso. Proceed with the regulation of the vertical axis (fig. C):  
- Extend the vertical splints and adapt the distance between the sternal and symphysis pelottes and to the length of the patient's torso;  
- Set the splint length in such a way that the horizontal section is at the height of the vertebra to which the lumbar pelotte is to be applied (roughly at the lumbar vertebra transition, or better, as recommended by the doctor);  
- Tighten the screws on the vertical axis and ensure that the orthosis is not uncomfortable when the patient is seated; if necessary adjust the vertical splint again;

3 Continue with the setting for the horizontal axis (fig. D):  
- Extend the horizontal splints to such an extent that the orthosis is adapted to the stomach circumference of the patient;  
- Tighten the screws and check the symmetry to the sagittal plane;

4 The orthosis can be adapted precisely to the shape of the patient's body:  
- Set the vertical bar (fig. E) so that the pressure exerted by the orthosis is directed to two points on the front of the body, the sternum and symphysis; if necessary re-position the sternal and symphysis pelottes (see point 2);  
- Set the horizontal fastenings (fig. F) in such a way that the pelvic bar rests on the hips of the patient.

5 Continue with the adjustment of the strap length:  
- Secure the hook to the round-headed screw (fig. H);  
- With the closure open, pull the strap through and secure with the Velcro® closure designed for this purpose;  
- If the strap is too long, it can be cut down accordingly and the Velcro® closure positioned in a y-shape (leave the strap as long as possible);  
- If necessary, pull the strap through again to obtain the desired lumbar pressure: optimum pressure is obtained when the hyperextension orthosis remains in position even with the closure opened (fig. G);  
- Secure the hyperextension orthosis with the help of the closure.

6 Ensure that the lumbar pelotte is positioned well against the spine (fig. I).

**PUTTING ON THE APPLIANCE FOR THE PATIENT**

1 Open the closure and release the strap from the round-headed screw.  
2 Place the hyperextension orthosis against the chest, pass the strap around the back of the body and secure to the round-headed screw (fig. L).  
3 Position the hyperextension orthosis in such a way that the sternal plate lies on the middle of the chest and the symphysis pelotte does not bother you when you are seated (fig. M).  
4 Secure the hyperextension orthosis using the closure (fig. N)

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