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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

## “dls” Modular spinal brace

### DECLARATION OF CONFORMITY

ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

### SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels.

Furthermore, we recommend wearing a garment to avoid direct contact with the skin.

If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product.

We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

### WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient.

To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

### SELECTION/SIZE

| Code              | REF.2701 |        |         |         |
|-------------------|----------|--------|---------|---------|
| Size              | S        | M      | L       | XL      |
| Underarm circ. cm | 75/85    | 85/95  | 95/105  | 105/115 |
| Waist circ. cm    | 65/80    | 80/95  | 95/110  | 110/125 |
| Pelvic circ. cm   | 70/85    | 85/100 | 100/115 | 115/130 |

### MAINTENANCE

- Washing instructions: wipe the plastic and synthetic leather parts with a soft sponge dipped in lukewarm water and mild soap. Dry with a cloth.
- Hand wash the product in tepid 30°C water with mild soap. Leave to dry away from direct heat.
- Periodically check if the screws need tightening.
- Consult an orthopaedic technician for the replacement of worn parts.

Do not dispose of the product or any of its components into the environment.

### INDICATIONS

- For serious dorsal-lumbar pathologies

### CONTRAINDICATIONS

Currently no known.

### FEATURES AND MATERIALS

- Combined structure
- Lightweight aluminum alloy frame
- Adjustable underarm bars on the upright rods
- Padded deltapectoral pads

### PUTTING ON THE APPLIANCE

- Adjust the upper part of the frame, shaping it as required.
- Insert the upper support into the rear pocket of the orthosis.
- Screw the axillary part, complete with pads and lining, to the upper part.
- Loosen the orthosis ties.
- Hook (for the women's version), or button (for the men's version) the orthosis at the front, starting from the bottom.
- Tighten the straps progressively starting from the bottom towards the top and put them in the corresponding loops. The adjustable buckles are self-blocking.

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