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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Arm and shoulder immobilizer with elbow support

DECLARATION OF CONFORMITY
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS
Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS
We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend not to wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS
The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician. The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements. The product must be fitted with the utmost care to ensure its efficacy, tolerability and correct function. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. The manufacturer is no longer liable if the product is used or adapted inappropriately. The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used.

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Item Size	REF.8409 universal
---------------------	------------------------------

Colour: grey - fits right and left

MANTENIMENTO

<input checked="" type="checkbox"/> Do not bleach	<input type="checkbox"/> No chemical cleaning
<input checked="" type="checkbox"/> Do not iron	<input type="checkbox"/> Do not tumble-dry
Washing instructions: hand wash at max. 30°C with neutral soap. Dry away from heat sources.	

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Stable humerus fractures
- Gleno-humeral luxation
- Distorsions of the shoulder
- Post surgery humerus fractures and gleno-humeral luxation treatment
- Post surgery after tenorrhaphy of long head of humeral biceps

CONTRAINDICATIONS
Currently no known

FEATURES AND MATERIALS

- Mixed cotton sling
- Immobilizer band made of soft padded fabric
- Shoulder straps with
- Velcro® fastening system
- Elbow pocket support

DONNING THE SUPPORT FOR THE FIRST TIME

- Attach the Velcro® fastenings to the ends of the two straps (A) as shown in the picture. If the straps are too long, you can cut off the excess part and then attach the Velcro® fastenings.
- Place the bent arm in the fabric sling (C).
- Pass the straps (B) over the shoulders.
- Cross over the straps at the back (D).
- Bring the ends of the straps (B) around to the front and secure them to the Velcro® loop fastening (E) sewn on the sling.
- Secure the swathe (F) to the Velcro® loop fastening (E) sewn on the sling.
- Wrap the swathe (F) around the torso and over the arm requiring immobilisation and then fasten it back on itself.

SUBSEQUENT WEAR

- Place the bent arm in the fabric sling (C).
- Pass the straps (B) over the shoulders.
- Cross over the straps at the back (D).
- Bring the ends of the straps (B) around to the front and secure them to the fastening.
- Wrap the swathe (F) around the torso and over the arm requiring immobilisation and then fasten it back on itself.

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