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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM**

Immobiliser for arm and shoulder with removable elbow pocket

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/EEC and 2007/47/EC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. Please read the composition of the product on the inside label carefully. You should not wear the device in the vicinity of an open flame.

WARNINGS


The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician. The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements. The product must be fitted with the utmost care to ensure its efficacy, tolerability and correct function. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. The manufacturer is no longer liable if the product is used or adapted inappropriately.

The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used.

SELECTION/SIZE

Code	REF: 8412
Size	universal size
Colour: black / Fits right and left	

MAINTENANCE

-  Do not bleach
-  No chemical cleaning
-  Do not iron
-  Do not tumble-dry
-  Washing instructions: hand wash at max. 30° C with neutral soap.
- Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Simple humeral fractures
- Glenohumeral luxation
- Distorsions of the shoulder
- Postoperative care of humeral fractures and interventions due to dislocations of the shoulder
- Postoperative care for tendon sutures of the long head of the biceps muscle

KONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Arm pocket made of padded fabric, hook and loop function on the outside, for optimum setting options
- Soft hook and loop shoulder straps
- Exclusive elbow support band made of mesh fabric
- Fixation bandage made of padded fabric, hook and loop function on the outside, for optimum setting options

PUTTING ON THE APPLIANCE

- Insert the arm of the shoulder to be immobilised into the arm pocket. Ensure that the orange label is outside and points upwards (fig. A).
- Place the straps over the shoulders (fig. B), cross them at the back (fig. C), then guide them to the front across the arm pocket and close them with the Velcro® closure (fig. D). The forearm has to be horizontal, parallel to the ground.
- Attach the fixation band using the corresponding Velcro® closure at the arm pocket (fig. E), guide around the chest (across the arm to be immobilised and under the opposite arm) and close the band again with the Velcro® closure (fig. F).
- Attach the elbow support band by fixing it on one side at the fixation band and on the other side at the arm pocket (fig. G).

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