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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM**

Arm and shoulder immobilizer with removable elbow pocket

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF: 8412
Size	universal
Colour: black Fits right and left	

MAINTENANCE

⊗ Do not bleach ⊗ No chemical cleaning

⚠ Do not iron ⚠ Do not tumble-dry

🧼 Washing instructions: hand wash at max. 30° C with neutral soap.

⚠ Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Stable humerus fractures
- Gleno-humeral luxation
- Distorsions of shoulder
- Post surgery humerus fractures and gleno-humeral luxation treatment
- Post surgery after tenorrhaphy of long head of humeral biceps

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Pocket sling made of
- Velcro padded fabric, perfectly adjustable
- Soft Velcro braces
- Exclusive mesh elbow-block band
- Fixation band made of padded fabric, perfectly adjustable

PUTTING ON THE APPLIANCE

- Insert the arm of the shoulder to be immobilised into the arm pocket. Ensure that the orange label is outside and points upwards (fig. A).
- Place the straps over the shoulders (fig. B), cross them at the back (fig. C), then guide them to the front across the arm pocket and close them with the Velcro® closure (fig. D). The forearm has to be horizontal, parallel to the ground.
- Attach the fixation band using the corresponding Velcro® closure at the arm pocket (fig. E), guide around the chest (across the arm to be immobilised and under the opposite arm) and close the band again with the Velcro® closure (fig. F).
- Attach the elbow support band by fixing it on one side at the fixation band and on the other side at the arm pocket (fig. G).

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