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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM

Brace for shoulder immobilization

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames. Do not apply in contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient.

Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation.

In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF.20413
Size	universal size
Colour	black
Fits right and left	

MAINTENANCE

☒ Do not bleach ☒ No chemical cleaning
☒ Do not iron ☒ Do not tumble-dry
☒ Washing instructions: hand wash at max. 30° C with neutral soap (a sponge is recommended); do not dry in the vicinity of heat sources (or, better, dry with a cloth)
Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Gleno-humeral luxation
- Distorsions of shoulder
- Post surgery

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Chest band made of velcro padded fabric, perfectly adjustable
- One arm band and two forearm bands made of Velcro® padded fabric to be freely positioned on chest band

APPLICATION

- Chest band: wrap around the chest at elbow height, and fasten with Velcro® (fig. 1a). If too long, cut off the excess part (fig. 1b).
- Arm band (orange trim): attach to the chest band using the Velcro® provided so that it is higher than the elbow (fig. 2a). Wrap it around the arm and fasten with Velcro® (fig. 2b).
- First forearm band (black trim): attach to the chest band using the Velcro® provided so that it is in the middle of the forearm when the elbow is bent (fig. 3a). Wrap it around the forearm and fasten with Velcro® (fig. 3b).
- Second forearm band (black trim): attach to the chest band using the Velcro® provided so that it is at hand height when the elbow is bent (fig. 4a). Wrap it around the hand and fasten with Velcro® (fig. 4b).

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Velcro® jest zastrzeżonym znakiem towarowym Velcro Industries B.V.