Axelimmobilisering ortos

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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM

Brace for shoulder immobilization

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames. Do not apply in contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation.

In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF.20414			
Size	S	M	L	XL
Forearm length* cm	28/34	34/42	42/48	48/58
Colour	black			

Fits right and left

*measurement taken from the elbow to the end of the clenched fist

MAINTENANCE

☒ Do not bleach ☒ No chemical cleaning
☒ Do not iron ☒ Do not tumble-dry
☒ Washing instructions: hand wash at max. 30° C with neutral soap (a sponge is recommended); do not dry in the vicinity of heat sources (or, better, dry with a cloth)
Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Stable humerus fractures
- Gleno-humeral luxation
- Distorsions of shoulder
- Post surgery humerus fractures and gleno-humeral luxation treatment
- Post surgery after tenorrhaphy of long head of humeral biceps

CHARACTERISTICS AND MATERIALS

- Made of cotton and polyester with elbow pocket support
- Chest immobilizer strap with central padded fabric and Velcro® fastening system
- Shoulder straps with Velcro® fastening system with central padded fabric

KONTRAINDICATIONS

Currently no known

APPLICATION

INITIAL FITTING

- 1 Undo the Velcro® fastenings on the chest strap and shoulder strap.
- 2 Place the bent arm in the fabric sling (fig. 1).
- 3 Pass the shoulder strap (black) over the opposite shoulder (fig. 2a), place the tip through the corresponding black ring, and fasten with Velcro® (fig. 2b).
- 4 Shorten the strap if necessary by cutting off the excess part (fig. 3a) and reattach the Y-Velcro® fasteners (fig. 3b). Check that the padded part remains centrally positioned to protect the neck.
- 5 Pass the (white) chest strap around the patient's chest, place the tip through the corresponding white ring (fig. 4a), and fasten with Velcro® (fig. 4b).
- 6 Shorten the strap if necessary by cutting off the excess part (fig. 5a) and reattach the Y-Velcro® fasteners (fig. 5b). Check that the padded part remains centrally positioned.
- 7 **SUBSEQUENT WEAR**
- 7 Undo the Velcro® fastenings on the chest strap and shoulder strap.
- 8 Place the bent arm in the fabric sling (fig. 10).
- 9 Pass the shoulder strap (black) over the opposite shoulder (fig. 11a), place the tip through the corresponding black ring, and fasten with Velcro® (fig. 11b).
- 10 Pass the (white) chest strap around the patient's chest, place the tip through the corresponding white ring (fig. 12a), and fasten with Velcro® (fig. 12b).

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