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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Shoulder brace with immobilizer in fixed abduction of 10°

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation.

The orthosis is made and to be used by only one patient.

Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF.90080	
Size	S	L
Forearm length cm*	< 38	> 39
Circumf. max. forearm cm	< 32	> 29
Colour	black	

specify right or left

*Measurement taken from elbow to end of clenched fist

MAINTENANCE

- Do not bleach
- No chemical cleaning
- Do not iron
- Do not tumble-dry
- Washing instructions: hand wash the cover (remove foam filling before doing so), the straps and cuff with warm water and a neutral soap. Rinse thoroughly. Do not dry in the vicinity of a heat source.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Post surgery (in abduction) to repair rotator cuff and soft tissue, glenoid labrum suture, capsular stabilization

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- New exclusive tensioning system:
 - Shoulder support (rotating buckles) on homolateral/injured shoulder to ensure its physiological posture, to prevent neckache.
 - Controlateral fixing strap with orange buckle for a rapid opening.
 - Back triangle straps for scapular stability
- Pad: high density foam with gauze fabric cover (to perfectly manage Velcro®) and fabric with silver fiber in contact with the body (anti-bacterial and odor-free)
- Length adjustable immobilizer
- Arm pocket support made of breathable material with Velcro® fastening system: easy elbow mobility.
- Adjustable abdominal strap with Velcro® fastening system and snap closing.
- Pad with front edge to improve stability. Tilted toward arm support for a better comfort
- Supporting strap for thumb and ball for rehabilitative purposes.

ADAPTING FOR THE DOCTOR / TECHNICIAN

- 1 Make sure that all straps have been pulled out fully and ONLY open the clip buckle of the abdominal strap (fig. A) and the orange hook closure (the front, contralateral strap) (fig. B). Also open the arm pocket and the strap used to immobilise the arm.
- 2 Guide the arm of the patient between the pad and the shoulder strap and position the pad at the side of the patient at elbow height on the side that is to be immobilised.
- 3 Close the clip buckle of the abdominal strap (fig. C) and tension using the appropriate Velcro® closure (fig. D). If the strap is too long: Remove the Y-shaped Velcro® closure end piece, shorten the strap to the appropriate length and reattach the Y-shaped Velcro® closure end piece.
- 4 Use the Velcro® closure (fig. E) to close the strap designed to immobilise the arm and to block the abduction movement of the shoulder. If it is too long, it can be shortened to the appropriate length and then reattached to the Velcro® closure on the rear of the pad.
- 5 Close the Velcro® closure of the arm pocket (fig. F) to secure the lower arm. If necessary, the size of the pocket can be adjusted by attaching the closure to a different position on the pad (using the Velcro® closure).
- 6 Position the thumb holder (with corresponding training ball) and secure to the pad with the Velcro® closure (fig. G).
- 7 Position the shoulder part as medially as possible (oriented towards the centre of the body) on the injured shoulder, insert the strap with the orange end piece under the opposite arm and attach using the hook closure on the pin of the shoulder part (fig. H). Secure the strap with the Velcro® closure to itself (fig. I), so that the position of the orange end piece corresponds to the marking on the strap in the same colour (this acts as a reference point when the patient is putting it on himself).
- 8 The contralateral straps are tensioned using the Y-shaped Velcro® closure (fig. K, L, M, N) in such a way that the shoulder part is secured and a triangle is formed on the rear which acts to stabilise the shoulder blade. If they are too long, shorten the straps and reattach the Y-shaped Velcro® end piece.
- 9 Position the padded cuff to protect the arm on the contralateral side (fig. O).
- 10 Tension the front vertical strap in such a way (fig. P) that the underarm is supported sufficiently. If the strap is too long: Remove the Y-shaped Velcro® closure end piece, shorten the strap to the appropriate length and reattach the Y-shaped Velcro® closure end piece.

REMOVING THE ORTHOSIS

- 1 Open the orange buckle (front contralateral strap), to do this unhook the buckle from the pin (if necessary, loosen the strap via the rubberised Velcro® closure end piece) (fig. A).
- 2 Open the arm pocket (fig. B) and the strap used to immobilise the arm (fig. C).
- 3 Open the clip buckle on the abdominal strap (fig. D) and remove the support carefully. Avoid sudden or extensive movements (fig. E).

PUTTING ON THE APPLIANCE FOR THE PATIENT

- 1 If not already open, open ONLY the clip buckle closure of the abdominal strap (fig. A) and the orange hook buckle (front contralateral strap), to do this unhook the buckle from the pin (fig. B). Also open the arm pocket and the strap used to immobilise the arm.
- 2 Insert the arm between the pad and the shoulder part and position the pad on the flank at elbow height (of the side to be immobilised).
- 3 Close the clip buckle of the abdominal strap (fig. C).
- 4 Use the Velcro® closure (fig. E) to close the strap designed to immobilise the arm and to block the abduction movement of the shoulder.
- 5 Close the Velcro® closure of the arm pocket (fig. F) to secure the lower arm.
- 6 Position the shoulder part as medially as possible (oriented towards the centre of the body) on the injured shoulder, insert the strap with the orange end piece under the opposite arm and attach using the hook closure on the pin of the shoulder part (fig. H). Secure the strap with the Velcro® closure to itself (fig. I), so that the position of the orange end piece corresponds to the marking on the strap in the same colour.

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Velcro® è un marchio di fabbrica registrato di Velcro Industries B.V.