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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Shoulder brace with immobilizer in fixed abduction

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully, do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient.

Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF. E110 10° - REF. E210 15°	REF. E111 10° - REF. E211 15°
Size	S M L	S L
Forearm length cm*	< 34 35/39 > 40	< 38 > 39

Colour: black - fits right and left
*Measurement taken from elbow to end of closed fist

Code	REF. E410
Size	universal

Colour: black - specify right or left

MAINTENANCE

-  Do not bleach
-  No chemical cleaning
-  Do not iron
-  No do tumble-dry
-  Washing instructions: Wash the cushion cover (after the foam inner padding has been removed), belts and the neck and underarm padding by hand with warm water (at max. 30°C) and neutral soap; rinse thoroughly. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Post surgery (in abduction) to repair rotator cuff and soft tissue, glenoid labrum suture, capsular stabilization

KONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Pad: high density foam with silver fiber in contact with the body
- Pad with front edge to improve stability; (only REF.410) Tilted toward arm support for a better comfort
- Protection neck sleeve
- External gauze nylon cover and sleeves with Velcro® fastening system
- Breathable arm pocket material
- Fixing straps with adjustable buckles with Velcro® fastening system
- Supporting strap for thumb

PUTTING ON THE APPLIANCE

ADAPTING FOR THE DOCTOR/TECHNICIAN

- 1 Undo the two side-release buckles. (REF. 410 ONLY) If necessary, detach the sling from the pillow and reposition it on the actual pillow so it matches up with the marking for the desired degree of abduction: orange=10°; white=15° (fig. A1). (EXCEPT FOR REF. 410) If necessary, detach the sling from the pillow and reposition it with Velcro® on the actual pillow on the same side as the shoulder requiring immobilisation. (fig. A2)
- 2 Position the pillow on the hip on the same side as the shoulder requiring immobilisation.
- 3 Fasten the waist strap with the side-release buckle (fig. B) and tighten with the corresponding Velcro® (fig. C). If it is too long, cut to size and reposition the Y-shaped Velcro® strap.
- 4 Carefully slide the forearm into the sling (fig. D1), fasten it with the Velcro® strap provided and (REF. 410 ONLY) adjust the length of the back strap (fig. D2).
- 5 Pass the vertical strap over the opposite shoulder and fasten the siderelease buckle (fig. E).
- 6 Tighten the shoulder strap with the corresponding Velcro®. If it is too long, cut to size and reposition the Y-shaped Velcro® strap.
- 7 Position the padded sleeve to protect the neck.
- 8 Position the thumb support (and exercise ball) and attach it to the pillow with the Velcro® (fig. F).

PUTTING ON THE APPLIANCE FOR THE PATIENT

Retain the configuration as adjusted by the doctor/technician. In order to remove the orthotic, open the arm fixation and unhook the push-in buckle fasteners. To put the orthotic on again, proceed as follows:

- 1 Position the pillow on the hip on the same side as the shoulder requiring immobilisation and rest the arm in the sling.
- 2 Fasten the waist strap with the side-release buckle (fig.B).
- 3 Close the sling with the Velcro fasteners and slide the thumb into its support.
- 4 Pass the vertical strap over the opposite shoulder and fasten the siderelease buckle (fig.E).

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