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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

## Shoulder brace with immobilizer in fixed abduction of 10°, underarm strap

### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

### SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

### WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

### SELECTION/SIZE

Code	REF.E112 10°	
Size	S	L
Forearm length cm*	< 38	>39

\*Measurement taken from elbow to end of closed fist

Colour: black - Right or left

### MAINTENANCE

- ⊗ Do not bleach
- ⊗ No chemical cleaning
- ⊗ Do not iron
- ⊗ Do not tumble-dry
- ⊗ Washing instructions: Wash the cushion cover (after the foam inner padding has been removed), belts and the neck and underarm padding by hand with warm water (at max. 30°C) and neutral soap; rinse thoroughly. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

### INDICATIONS

- Post surgery (in abduction) to repair rotator cuff and soft tissue, glenoid labrum suture, capsular stabilization

### CONTRAINDICATIONS

Currently no known

### CHARACTERISTICS AND MATERIALS

- High density foam rubber with gauze fabric cover (to perfectly manage Velcro®) and fabric with silver fiber in contact with the body
- Pad with front edge for an excellent stability
- External gauze nylon cover and sleeves with Velcro® fastening system
- Breathable arm pocket material
- Sleeve for neck protection
- Underarm strap to avoid neck pressure and increase rotatory stability
- Fixing straps with adjustable buckles with Velcro® fastening system
- Ball for rehabilitative purposes.
- Supporting strap for thumb
- Abduction of 10°

### PUTTING ON THE APPLIANCE

#### ADAPTING FOR THE DOCTOR/TECHNICIAN

- 1 If necessary, remove the arm sling from the cushion and attach to the cushion using the Velcro® fastener so that the orthotic lies to the left or the right-hand side.
- 2 Place the cushion laterally at the side to which the shoulder should be fixed, close the waist belt with the push-in buckle fastener (fig. A) and tighten with the corresponding Velcro® fastener (fig. A). Should the belt be too long, then it can be shortened. To do this, cut the belt to the appropriate length and re-attach the Y-shaped Velcro® fastener end.
- 3 Carefully insert the forearm into the arm fixation (fig. B); then close the arm fixation with the respective Velcro® fasteners at the wrist (fig. C) and the forearm (fig. D).
- 4 Insert the arm on the opposite side from the shoulder to be fixed through the loop created by the shoulder and underarm belts (fig. E) and close the orthotic with the push-in buckle fastener on the arm fixation (fig. F).
- 5 Close all the belts pertaining to the healthy shoulder with the respective Velcro® fasteners (fig. G); Back belt, underarm belt, shoulder belt and front belt. If the belt is too long, it can be shortened; to do this, cut the belt to the appropriate length and re-attach the Y-shaped Velcro® fastener end.
- 6 Check the placement of the underarm and neck padding; re-position the padding if required (fig. H).
- 7 Attach the thumb support (and the softball for rehabilitation) onto the cushion with the Velcro® fastener (fig. I).

#### PUTTING ON THE APPLIANCE FOR THE PATIENT

Retain the configuration as adjusted by the doctor/technician. In order to remove the orthotic, open the arm fixation and unhook the push-in buckle fasteners. To put the orthotic on again, proceed as follows:

- 1 Place the cushion laterally at the side on which the shoulder is to be fixed; then place the arm into the breathable arm fixation.
- 2 Close the waist belt with the push-in buckle fastener (fig. L).
- 3 Close the arm fixation with the respective Velcro® fasteners and insert the thumb into the thumb support.
- 4 Insert the arm on the opposite side from the shoulder to be fixed through the loop created by the shoulder and underarm belts (fig. E); Close the orthotic with the push-in buckle fastener on the arm fixation (fig. F).

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