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**PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

## Shoulder brace with immobilizer in external rotation of 15°/20°

### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

### SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

### WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

### SELECTION/SIZE

Code	REF.84502
Size	universal size
Colour:	black / Fits right and left
Option:	REF.84500/A - Hand and wrist support

### MAINTENANCE

- Do not bleach  No chemical cleaning
- Do not iron  Do not tumble-dry
- Washing instructions: wash the cover (after having removed the foam padding), the belt and the cuffs by hand in lukewarm water with neutral soap, rinse thoroughly. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

### INDICATIONS

- Post-operative treatment of rotator cuff reconstructions
- Nerve injuries

### CONTRAINDICATIONS

Currently no known

### CHARACTERISTICS AND MATERIALS

- High density foam pad
- Cover, chest straps and sleeves made of paired fabric with external gauze nylon to be fastened with Velcro®
- Shoulder strap with Velcro® fastening system
- Abduction of 15° or 20°

### APPLICATION

#### FIRST APPLICATION BY THE DOCTOR/TECHNICIAN

- 1 Position the pillow in order to ensure the desired abduction (15° or 20°). Place the pillow on the side, above the iliac crest. Tighten the band around the waist and attach it to the pillow with the Velcro® (fig. A).
- 2 If the rod is being used (optional REF. 84500/A), attach it to the pillow with the hand support placed in front.
- 3 Position and attach the humerus stabilizing cuff to the pillow using the Velcro®, followed by the wrist cuff.
- 4 Close the humerus stabilizing cuff (fig. B), followed by the wrist cuff (fig. C).
- 5 Should the position of the arm be unnatural, reposition the cuffs on the pillow.
- 6 Attach the shoulder strap to the front (fig. D) and back (fig. E) of the pillow.

### SUBSEQUENT APPLICATIONS

- 1 Place the pillow on the side above the iliac crest.
- 2 Tighten the band around the waist and attach it to the pillow with the Velcro® (fig. A).
- 3 Close the humerus stabilizing cuff (fig. B), followed by the wrist cuff (fig. C).
- 4 Attach the shoulder strap to the front (fig. D) and back (fig. E) of the pillow.

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