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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Shoulder brace with immobilizer in abduction of 45°/70°

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation.

The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation.

In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF. 84500
Size	universal
Colour	black

Fits right and left
Option: REF. 84500/A - Hand and wrist support

MAINTENANCE

⊗ Do not bleach ⊗ No chemical cleaning

⊗ Do not iron ⊗ Do not tumble-dry

⊗ Washing instructions: wash the cover (after having removed the foam padding), the belt and the cuffs by hand in lukewarm water with neutral soap, rinse thoroughly. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Tuberosity fractures with slight dislocation
- Adhesive capsulitis after mobilization
- Post-surgery of endo and arthroplasty of shoulder
- Post-surgery of rotator cuff injury
- Nerve injury

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- High density foam pad
- Cover, chest strap and sleeves made of paired fabric with external gauze nylon to be fastened with Velcro®
- Shoulder strap with Velcro® fastening system
- Abduction of 45° or 70° is possible
- Protection neck sleeve

ADAPTING FOR THE DOCTOR / TECHNICIAN

- 1 Position the cushion in such a way that the desired abduction (45° or 70°) is achieved. Place the cushion sideways over the iliac crest. Put the belt around the waist and fasten with the Velcro® closure (fig. A).
- 2 If the hand and wrist support is used (optional REF. 84500/A), it has to be attached to the cushion in such a way that the hand support points forward (fig. B).
- 3 Position the upper arm cuff (fig. C) and subsequently the wrist cuff (fig. D) and fasten with the Velcro® closure.
- 4 Close the upper arm cuff and then the wrist cuff. In case the shoulder position is unnatural, change the position of the cuffs on the cushion.
- 5 Fix the shoulder strap on the rear side (fig. E) and on the front side of the cushion (fig. F), guide through the appropriate ratchet buckle and tension. Ensure that the neck cushion is positioned correctly.

PUTTING ON THE APPLIANCE FOR THE PATIENT

- 1 Place the cushion sideways over the iliac crest.
- 2 Pull the belt around the waist and fix it at the cushion using the Velcro® closure (fig. A).
- 3 Close the upper arm cuff (fig. C) and subsequently the wrist cuff (fig. D).
- 4 Attach the shoulder strap on the rear side (fig. E) and on the front side of the cushion (fig. F), guide it through the appropriate ratchet buckle and tension. Ensure that the neck cushion is positioned correctly.

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