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Abduction of 45°/60°  
(depending on inflation level)

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#### SCOPE OF DELIVERY

- 1 Adhesive Velcro® straps with closures and three-cornered eyelets
- 2 Inflatable inner cushion (optionally available with pump REF. 84505/P)
- 3 Cover with buckles and Velcro® closure
- 4 Wedge-shaped arm support and handsupport
- 5 Abdominal strap
- 6 Lower arm securing straps

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REF. 84503

Shoulder abduction wedge

DESIGNED AND  
DEVELOPED IN  
SWITZERLAND

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
Abduction of 45°/60° (depending on inflation level)

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SCOPE OF DELIVERY

- 1 Adhesive Velcro® straps with closures and three-cornered eyelets
2 Inflatable inner cushion (optionally available with pump REF. 84505/P)
3 Cover with buckles and Velcro® closure
4 Wedge-shaped arm support and handsupport
5 Abdominal strap
6 Lower arm securing straps

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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM
Shoulder abduction wedge

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR).

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient.

SELECTION/SIZE

Table with columns: Code, Size, Colour and corresponding values: REF. 84503, Universal, black.

fits right and left
OPTION: REF. 84505/P with pump

MAINTENANCE

- ☒ Do not bleach
☒ No chemical cleaning
☒ Do not iron
☒ Do not tumble-dry
Washing instructions: wash the cover (after having removed the foam padding), the belt and the cuffs by hand in lukewarm water with neutral soap, rinse thoroughly.

INDICATIONS

- Tuberculum-majus fracture with minor dislocation
• Treatment of capsulitis adhaesiva (frozen shoulder) after mobilisation
• Post-operative care of endoprosthesis and arthroprosthesis of the shoulder
• Post-operative care after surgery to the rotator cuff
• Nerve damages

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Structure of the cushion: Cylinder-shaped, inflatable inner cushion made from robust plastic
• Exterior made using gauze fabric with adhesive surfaces for the Velcro® closure
• Shoulder and pelvic strap made from nylon gauze with Velcro® closure
• Pad with sandwich structure with interior padding
• Wedge-shaped axilla support made from coated foam for improved stability and positioning of the cushion
• Abduction can be regulated by modifying the inflation of the cushion
• Lower arm and hand support

ADAPTING FOR THE DOCTOR / TECHNICIAN

- 1 Insert the inner cushion into the cover which is usable with Velcro® (if the adhesive surface is already attached, remove it). Inflate the inner cushion to the required volume. The cushion must be positioned against the side of the patient to better ensure stability.
2 Insert the arm of the healthy side through the ring formed by the shoulder and axilla strap so that the eyelet that is affixed to the adjustable strap faces forwards after the arm has been moved to the correct position (fig. A).
3 The eyelets sewn to the straps must, on the other hand, be on the back (fig. B).
4 Secure the cushion with the appropriate buckles to the straps (fig. C).
5 Position the hand support in the most convenient holding position and secure it to the cover (fig. D).
6 Grasp the hand support and secure the two lower arm straps (fig. E) in accordance with the doctor's instructions.
7 Affix the wedge-shaped arm support under the medial (inner) side of the arm. (There should be no body contact, if possible.)

PUTTING ON THE APPLIANCE FOR THE PATIENT

- 1 Insert the arm of the healthy shoulder through the ring formed by the shoulder and axilla strap so that the eyelet that is affixed to the adjustable strap faces forwards after the arm has been moved to the correct position (fig. A).
2 The eyelet sewn to the straps must, on the other hand, be on the back (fig. B).
3 Secure the cushion with the appropriate buckles to the straps (fig. C).
4 Grasp the hand support and secure the two lower arm straps (fig. D).
5 Affix the wedge-shaped arm support under the medial (inner) side of the arm. (There should be no body contact, if possible.)