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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM
Shoulder abduction wedge

DECLARATION OF CONFORMITY
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS
Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS
We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS
It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

Code	REF. 84503
Size	Universal
Colour	black

fits right and left
OPTION: REF. 84505/P with pump

MAINTENANCE
☒ Do not bleach ☒ No chemical cleaning
☒ Do not iron ☒ Do not tumble-dry
☒ Washing instructions: wash the cover (after having removed the foam padding), the belt and the cuffs by hand in lukewarm water with neutral soap, rinse thoroughly. Dry away from heat sources. Do not dispose of the product or any of its components into the environment.

INDICATIONS
• Tuberculum-majus fracture with minor dislocation
• Treatment of capsulitis adhaesiva (frozen shoulder) after mobilisation
• Post-operative care of endoprosthesis and arthroprosthesis of the shoulder
• Post-operative care after surgery to the rotator cuff
• Nerve damages

CONTRAINDICATIONS
Currently no known

CHARACTERISTICS AND MATERIALS
• Structure of the cushion: Cylinder-shaped, inflatable inner cushion made from robust plastic
• Exterior made using gauze fabric with adhesive surfaces for the Velcro® closure
• Shoulder and pelvic strap made from nylon gauze with Velcro® closure
• Pad with sandwich structure with interior padding
• Wedge-shaped axilla support made from coated foam for improved stability and positioning of the cushion
• Abduction can be regulated by modifying the inflation of the cushion
• Lower arm and hand support

ADAPTING FOR THE DOCTOR / TECHNICIAN
1 Insert the inner cushion into the cover which is usable with Velcro® (if the adhesive surface is already attached, remove it). Inflate the inner cushion to the required volume. The cushion must be positioned against the side of the patient to better ensure stability. Use the appropriate adhesive surface to close the cover (and conceal the valve).
2 Insert the arm of the healthy side through the ring formed by the shoulder and axilla strap so that the eyelet that is affixed to the adjustable strap faces forwards after the arm has been moved to the correct position (fig. A).
3 The eyelets sewn to the straps must, on the other hand, be on the back (fig. B).
4 Secure the cushion with the appropriate buckles to the straps (fig. C). Regulate the length of the straps so that the arm lies on the cushion with the desired angle of abduction of the shoulder.
5 Position the hand support in the most convenient holding position and secure it to the cover (fig. D). Provisionally attach both straps for the lower arm.
6 Grasp the hand support and secure the two lower arm straps (fig. E) in accordance with the doctor's instructions Note: Ensure that the head of the humerus is not pushed upwards.
7 Affix the wedge-shaped arm support under the medial (inner) side of the arm. (There should be no body contact, if possible.) Secure the wedge using the appropriate strap (fig. F).

PUTTING ON THE APPLIANCE FOR THE PATIENT
1 Insert the arm of the healthy shoulder through the ring formed by the shoulder and axilla strap so that the eyelet that is affixed to the adjustable strap faces forwards after the arm has been moved to the correct position (fig. A).
2 The eyelet sewn to the straps must, on the other hand, be on the back (fig. B).
3 Secure the cushion with the appropriate buckles to the straps (fig. C).
4 Grasp the hand support and secure the two lower arm straps (fig. D).
Note: Ensure that the head of the humerus is not pushed upwards.
5 Affix the wedge-shaped arm support under the medial (inner) side of the arm. (There should be no body contact, if possible.) Secure the wedge using the appropriate strap (fig. E).