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
 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Shoulder brace with immobilizer in external rotation of 15° or 30°

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care.

Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item right	REF. E113 15°	REF. E115 30°		
Item left	REF. E114 15°	REF. E116 30°		
Size	S	M	L	
Forearm length cm*	< 34	35/39	> 40	
Colour	black			

*Measurement taken from elbow to end of closed fist
Specify right or left

MAINTENANCE

 Do not bleach  No chemical cleaning

 Do not iron  Do not tumble-dry

 Washing instructions: Wash the cushion cover (after the foam inner padding has been removed), belts and the neck and underarm padding by hand with warm water (at max. 30°C) and neutral soap; rinse thoroughly. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Post surgery (external rotation) to repair rotator cuff and soft tissue, glenoid labrum suture, capsular stabilization
- Conservative treatment of first glenohumeral dislocations

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- High density foam pad with gauze fabric cover made of silver fiber in contact with the body
- Cover and sleeves made of external gauze nylon to be fastened with Velcro®
- Adjustable straps with buckles and Velcro® fastening system
- Protection neck sleeve
- Possible external rotation of 15° or 30°
- Supporting strap for thumb
- Ball for rehabilitative purposes.
- Arm pocket support made of breathable material

PUTTING ON THE APPLIANCE

ADAPTATION FOR THE DOCTOR/TECHNICIAN

- 1 Position the cushion on the side corresponding to the shoulder to be immobilised and place the arm in the appropriate breathing support pocket (fig. A). The arm-support pocket for the stabilisation of the limb may be repositioned on the cushion (using the Velcro®), if necessary.
- 2 Close the waistband with the large click-fasten buckle (fig. B) and tighten it with the corresponding Velcro® (fig. C). If it is too long, it may be shortened by cutting it to size and repositioning the Velcro® in a Y shape.
- 3 Close the arm-support pocket with the appropriate Velcro® strips (fig. D).
- 4 Position the thumb support (and the associated rehabilitation ball), using the Velcro® to attach it to the cushion (fig. E).
- 5 For correct extra-rotation positioning, position yourself in front of a table and make sure that the white line on the upper part of the cushion is parallel to the table.
- 6 Slip the opposite arm into the triangle created by the armpit and the clavicle strap, so that the white marked strap is above the arm and the white-marked one passes below the armpit (fig. F).
- 7 Close the guide by attaching the small click-fasten buckle (fig. G).
- 8 Tighten the (orange) shoulder strap with the appropriate Velcro® (fig. H). If it is too long, it may be shortened by cutting it to size and repositioning the Velcro® in a Y shape.
- 9 Tighten the (white) under-arm strap with the appropriate Velcro® (fig. I). If it is too long, it may be shortened by cutting it to size and repositioning the Velcro® in a Y shape.
- 10 Check and if necessary reposition the sleeve padding at the level of the neck and under the armpit.
- 11 Tighten the rear strap with the appropriate Velcro® (fig. K). If it is too long, it may be shortened by cutting it to size and repositioning the Velcro® in a Y shape.

PUTTING ON THE APPLIANCE FOR THE PATIENT

- Maintain the configuration set when using for the first time. The appliance may be removed by opening the arm-support pocket and unfastening the closure buckles. The support is reapplied as follows:
- 1 Position the cushion on the side corresponding to the shoulder to be immobilised and place the arm in the appropriate respiring support pocket (fig. A).
 - 2 Close the waistband with the large click-fasten buckle (fig. B).
 - 3 Close the arm-support pocket with the appropriate Velcro® strips (fig. D) and slip the thumb into the appropriate support.
 - 4 For correct extra-rotation positioning, position yourself in front of a table and make sure that the white line on the upper part of the cushion is parallel to the table.
 - 5 Slip the opposite arm into the triangle created by the armpit and the clavicle strap, so that the orange marked strap is above the arm and the white-marked one passes below the armpit (fig. F).
 - 6 Close the guide by attaching the small click-fasten buckle (fig. G).

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Velcro® è un marchio di fabbrica registrato di Velcro Industries B.V.