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REF-20726 - Functional shoulder brace

DECLARATION OF CONFORMITY
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). The instructions have been drawn up in application of the aforementioned Regulation. They are prepared in accordance with the applicable standards and ensure safe use of the medical device.

TRADEMARKS OF MATERIALS
Velcro® is a registered trademark of Velcro Industries B.V.
PRECAUTIONS
We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. In doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in contact with open wounds.

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Specify right or left

MAINTENANCE
 Do not bleach No chemical cleaning
 Do not iron Do not tumble-dry
 Washing instructions: Hand wash at max. 30° C with neutral soap (a sponge is recommended); do not dry in the vicinity of heat sources.
Do not dispose of the product or any of its components into the environment.

- INDICATIONS**
- Conservative treatment of instability following anterior or posterior dislocation of the glenohumeral joint
 - protection of the shoulder during at-risk sports or occupational activities while waiting for surgery or functional stabilisation.
 - Use for shoulder rehabilitation treatment of painful hyperlaxity or functional instability
 - Use three months after surgery to repair rotator cuff when starting work or sport, including when managing massive ruptures of rotator cuff
 - Helpful in managing the rehabilitation of scapular dyskinesia.

CONTRAINDICATIONS
Currently no known

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- PUTTING ON THE APPLIANCE - INITIAL FITTING**
- Undo the buckle and undo the Velcro® fastening on the arm.
 - Slide the arm into the tubular part, keeping the buckle towards the top. Pull the brace up the arm until it almost touches the neck and the shoulder (fig. 1).
 - C Attach the arm closure to Velcro® so that it fits well but without blocking circulation. (fig.2)
 - D Hold the rear fastening strap, pass it under the opposite armpit (fig.3a) and fasten the buckle at the front (fig.3b).
 - E Adjust the length of the fastening strap at the front (fig.4a) and at the back (fig.4b), so that the padding is centred under the armpit: cut off the excess part and secure the strap with the figure (fig.4b) fastener.
 - F Tighten and fasten the two elasticated straps as prescribed by a doctor (both straps can be passed through the upper loop or one through the lower loop and one through the upper loop, depending on the patient's needs).
 - G Move the arm a little a few times to check that the brace is very stable. If necessary, adjust the underarm strap again.

- SUBSEQUENT WEAR**
- Undo the buckle and turn the end part of the arm socket so that the non-slip material is external and does not make it difficult to insert the brace.
 - Slide the arm into the tubular part, keeping the buckle towards the top and the elasticated arm closure towards the bottom. Pull the brace up the arm until it almost touches the neck and the shoulder is completely covered. (fig.10)
 - J Hold the rear fastening strap, pass it under the opposite armpit (fig.11a) and fasten the buckle at the front (fig.11b).
 - K Tighten and fasten the two elasticated straps:
 - fasten the upper strap in the lower part of the arm socket, as indicated in the figure (fig.12a)
 - Crossover and fasten the lower strap (marked in orange) to the upper Velcro® (also marked in orange), situated near the buckle (fig.12b)

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