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EN FR IT PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

**REF.725 Humeral fracture brace**  
**REF.726 Shoulder cap humeral brace**

**DECLARATION OF CONFORMITY**

ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

**TRADEMARKS OF MATERIALS**

Velcro® is a registered trademark of Velcro Industries B.V.

**SAFETY PRECAUTIONS**

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

**WARNINGS**

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

**SELECTION/SIZE**

Code	REF.725 - REF.726				
	XS	S	M	L	XL
Größe					
REF.725 Proximal circ. humerus cm	-	23/28	28/33	33/38	38/43
REF.726 Proximal circ. humerus cm	20/25	23/28	28/33	33/38	38/43
REF.725 Distal circ. humerus cm	-	18/23	23/28	28/33	33/38
REF.726 Distal circ. humerus cm	15/20	18/23	23/28	28/33	33/38
REF.725 Brace length cm	-	19	20	21	23
REF.726 Brace length* cm	34	34	34	34	34
Colour	white				

Fits right and left  
 \* you may shorten by cutting it

**MAINTENANCE**

- ☒ Do not bleach ☒ No chemical cleaning
- ☒ Do not iron ☒ Do not tumble-dry
- ☒ Instructions for washing: wipe with a sponge dipped in lukewarm water and mild soap. Dry with a cloth. Consult an orthopaedic technician for the replacement of worn parts. Do not dispose of the product or any of its components into the environment.

**INDICATIONS**

- Stable fractures of medial distal third of humerus
- Post-surgery brace after humerus fractures
- Fractures due to non-operative bone metastasis (Only REF. 726)

**CONTRAINDICATIONS**

Currently no known

**CHARACTERISTICS AND MATERIALS**

- Structure made of Polyethylene (Only REF. 725)
- Structure made of Polyethylene (Only REF. 726)
- Internal pad made of close cells against moisture
- Extension cap to shoulder-humeral joint for better control of movements (Only REF. 726)
- Velcro® straps
- Ventilation holes
- Two tubular sleeves included

**PUTTING ON THE APPLIANCE REF.725**

- 1 Slip the protective sock onto the upper arm.
- 2 Don the brace and tighten the Velcro® straps. The open side must correspond with the inner side of the arm.
- 3 Don the arm support belt.

**PUTTING ON THE APPLIANCE REF.726**

- 1 If necessary, mould the plastic parts with the help of a heat gun.
- 2 Slip the protective sock onto the upper arm.
- 3 Don the brace and tighten the Velcro® straps.
- 4 Fasten the shoulder cap securely with the Velcro® chest strap.

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