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BREATHABLE

 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Standard arm sling

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician.

The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements. The product must be fitted with the utmost care to ensure its efficacy, tolerability and correct function.

Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician.

Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. The manufacturer is no longer liable if the product is used or adapted inappropriately.

The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals.


In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used.

SELECTION/SIZE

Item	REF.8404			
Size	S	M	L	XL
Pocket length cm	34	42	48	58
Colour	black			

fits right and left

MAINTENANCE

-  Do not bleach
-  No chemical cleaning
-  Do not iron
-  Do not tumble-dry

 Washing instructions: hand wash at max. 30° C with neutral soap rinse thoroughly.

Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Forearm and wrist support
- Support when cast
- Use in case of acromioclavicular (after 1° subluxation) or glenohumeral rest is required
- Use in case of upper arm immobilization with elbow in flexed position

CONTRAINDICATIONS

Currently no known

FEATURES AND MATERIALS

- Structure made of cotton blend
- Padded MTP Shoulder strap with velveteen protection
- Velcro® fastening system

PUTTING ON THE APPLIANCE

- 1 Place the bent arm in the fabric sling (fig.A)
- 2 Place the thumb through the thumbhole if necessary. (fig.B)
- 3 Pass the shoulder strap over the opposite shoulder to the arm which requires support (fig.C)
- 4 Pass the strap through the corresponding ring (fig. D) and pull as tightly as required (fig.E). The strap may be shortened using the Velcro® at the back

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