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
PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

### Medial epicondylitis strap providing compression

#### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR) herewith. The following instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

#### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

#### SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body which present ferite, gonfiori or tumefactions. It is advisable to not over tighten the device so as to avoid creating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels.

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#### AUSSAHL/GRÖSSEN

Code	REF.91020	XS	S	M	L	XL
Size		20/22	23/25	26/28	29/31	32/34
Max forearm circumference cm		black				

Fits right and left

#### PFLEGE

- Nicht bleichen
- No chemical cleaning
- Nicht bügeln
- Do not tumble-dry
- Washing instructions: **Pölsterung:** Hand wash in lukewarm water (30°C) with mild soap; rinse thoroughly. Leave to dry away from direct heat.
- Support strap:** wipe with a sponge dipped in warm water (30°C) with mild soap.

Do not dispose of the product or any of its components into the environment.

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#### CONTRAINDICATIONS

There are no known contraindications at this time.

#### FEATURES AND MATERIALS

- Rigid polypropylene C-shape design
- Velcro® fastening system
- Soft pad with double-lined MTP
- Pelote for targeted compression

#### PUTTING ON THE APPLIANCE

- 1 (Initial fitting only) Attach the pressure pad with Velcro® to the widest part of the padding; the long part of the pressure pad should be positioned at right angles to the long side of the strap.
- 2 Position the strap with the C-shaped part on the side so that the pressure pad is three or four fingers' width below the epicondyle (fig. A), and align it with the extensor muscles of the hand.
- 3 Tighten the strap and fasten it with the Velcro® (fig. B). Check that the pressure pad is in the correct position.
- 4 If false, the brace is too tight or causes a tingling sensation, loosen to relieve the pressure by adjusting the fastening.

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