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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

## Medial epicondylitis strap providing compression

### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

### SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/ or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

### WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

### SELECTION/SIZE

Code	REF.91020					
Size	XS	S	M	L	XL	
Max forearm circumference cm	20/22	23/25	26/28	29/31	32/34	
Colour	black					

Fits right and left

### MAINTENANCE

Do not bleach  No chemical cleaning

Do not iron  Do not tumble-dry

Washing instructions: **Padding:** Hand wash in lukewarm water (30°C) with mild soap;

rinse thoroughly. Leave to dry away from direct heat.

**Support strap:** wipe with a sponge dipped in warm water (30°C) with mild soap.

Dry with a cloth.

Do not dispose of the product or any of its components into the environment.

### INDICATIONS

- Severe and chronic epicondylitis and medial epicondylitis
- Muscle tensions
- Preventing inflammation during sports

### CONTRAINDICATIONS

There are no known contraindications at this time.

### FEATURES AND MATERIALS

- Rigid polypropylene C-shape design
- Velcro® fastening system
- Soft pad with double-lined MTP
- Targeted compression

### PUTTING ON THE APPLIANCE

- 1 (Initial fitting only) Attach the pressure pad with Velcro® to the widest part of the padding; the long part of the pressure pad should be positioned at right angles to the long side of the strap.
- 2 Position the strap with the C-shaped part on the side so that the pressure pad is three or four fingers' width below the epicondyle (fig. A), and align it with the extensor muscles of the hand.
- 3 Tighten the strap and fasten it with the Velcro® (fig. B). Check that the pressure pad is in the correct position.
- 4 If it feels too tight or causes a tingling sensation, loosen to relieve the pressure by adjusting the fastening.

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Velcro® jest zastrzeżonym znakiem towarowym Velcro Industries B.V.