# episystem

## REF. P1022

Multifunktionsbandage mit 2 abnehmbaren Druckpelotten  
Multi-functional bandage with 2 removable pressure pads  
Bandage multifonctions avec 2 pelotes de pression amovibles  
Многофункциональный ремешок с двумя типами съемных прижимов  
Cinturino multifunzionale dotato di 2 tipi di pressori rimovibili

## DEUTSCH

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## ENGLISH

Please read the instructions carefully and thoroughly

### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/EEC and 2007/47/EC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

### SELECTION/SIZE

Universal size · fits right and left

Mark no.	1	2	3	4
Circumference Forearm (cm)	20	25	30	35

A= Reference point for the alignment of the marks

### MAINTENANCE

- Do not bleach
- No chemical cleaning
- Do not iron
- Do not tumble-dry

Washing instructions: Remove the pads. Hand wash the bandage, the cover of the cooling pad and the pressure pad at max. 30°C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro® Industries B.V.

### PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. The device should not be worn in the vicinity of open flames.

### WARNINGS

The product may only be used in the event of medical indication or when prescribed by a physician and must be applied by an orthopaedic technician, the competent authority figure for the application as well as for information on the safe use, according to individual needs. To ensure the effectiveness, tolerability, and correct operation, it is necessary that the device is put on with the utmost care. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. Use of the product is recommended for a single patient. With hypersensitive patients, direct contact with skin may cause redness or irritation. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician.

### WARRANTY

This warranty is governed by the statutory provisions of the country in which the product was purchased. Where relevant that country's regulations concerning warranties and the relationship between the seller and buyer shall apply. If you believe that the terms of the warranty may be applicable, first contact the retailer where the product was purchased. Orthoservice will repair or replace all or part of the device and its accessories if defects in materials and workmanship arise during the period of six months from the date of purchase.

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## ENGLISH

- FEATURES AND MATERIALS**
- Rigid polypropylene C-shape design
  - Elasticated strap and Velcro® fastening
  - Soft double-lined MTP padding
  - Pressure pad for targeted pressure

- INDICATIONS**
- Acute or chronic epicondylitis and epitrochleitis of the elbow
  - Muscle contractions
  - Prevention of inflammatory processes during sporting activities
- CONTRAINDICATIONS**  
There are no known contraindications at this time.

**ADAPTING FOR THE DOCTOR/TECHNICIAN**

- 1 Insert the pressure pad or the cold pressure pad in the area provided (fig. A/B).
- 2 Widen the bandage (without pulling it out from the ring) and fit it (pull up from the hand towards the arm) (fig. C).
- 3 Place the bandage 3 or 4 fingers below the epicondyle (fig. D) and turn until the pressure pad is located on the extensor muscles of the hand (fig. E).
- 4 Open the Velcro® closure and tighten the bandage until the white line on the inside of the ring appears (connected to the elastic band (fig. F).
- 5 If required, remove the bandage to adjust its length: (fig. G/H)
  - a). Remove the Y-shaped Velcro® closure
  - b). Cut the protruding part
  - c). Re-attach the Y-shaped Velcro® closure.
- 6 Put on the bandage again and repeat steps 2 to 4, then close the bandage with the Velcro® closure (fig. I).
- 7 Check that the pressure pad is in the correct position (point 3).
- 8 Should the bandage disturb too much when wearing it or cause a tingling sensation, slightly loosen the closure.

**PUTTING ON THE APPLIANCE FOR THE PATIENT**

- 1 Widen the bandage (without pulling it out from the ring) and fit it (pull up from the hand towards the arm) (fig. C).
- 2 Place the bandage 3 or 4 fingers below the epicondyle (fig. D) and turn until the pressure pad is located on the extensor muscles of the hand (fig. E).
- 3 Open the Velcro® closure and tighten the bandage until the white line on the inside of the ring appears (connected to the elastic band) (fig. F).
- 4 Close the bandage with the Velcro® closure (fig. I).
- 5 Check that the pressure pad is in the correct position (point 2).
- 6 Should the bandage disturb too much when wearing it or cause a tingling sensation, slightly loosen the closure.

**COOLING**

- 7 Take the gel cushion out of its velour cover.
- 8 Place the cushion into the freezer at least 2-3 hours prior to its application (NOTE: Do NOT put the velour cover in the freezer).
- 9 Shortly before applying the gel cushion insert it in the velour cover (fig. L).
- 10 Replace the pressure pad that was used during sports with the gel cushion (fig. M).
- 11 Attach the bandage in the same way as for sporting activities. Ensure that the gel cushion is exactly in the same position as the first pressure pad.
- 12 Wear the bandage for the duration prescribed by the physician.

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