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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM

Epicondylitis strap with pressure pad

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/ EEC and 2007/47/EC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.




PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. If in doubt about how to apply, contact an orthopaedic technician. The device should not be worn in the vicinity of open flames. People prone to dermatitis or allergies to products resulting from the neoprene must not use this brace.

WARNINGS

The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician. The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements. The product must be fitted with the utmost care to ensure its efficacy, tolerability and correct function. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. The manufacturer is no longer liable if the product is used or adapted inappropriately. The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used.

MAINTENANCE

-  Do not bleach
-  No chemical cleaning
-  Do not iron
-  Do not tumble-dry
-  Washing instructions: Remove the pads. Hand wash the bandage, the cover of the cooling pad and the pressure pad at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

SELECTION/SIZE

Code	REF. 4320	
Size	S	L
Circumference Forearm (cm)	< 27	> 27
Colour	beige	

fits right and left

FEATURES AND MATERIALS

- Double-lined neoprene with leatherette support
- Velcro® fastening with a nylon loop
- Nonslip fabric pressure pad

INDICATIONS

- Acute or chronic tennis elbow and golfer's elbow
- Muscle contractions
- Prevention of inflammation during sports activities

CONTRAINDICATIONS

There are no known contraindications at this time.

PUTTING ON THE APPLIANCE

- 1 Attach the pressure pad on the bandage in the area Velcro® provided; the pressure pad has to be placed with the long part perpendicular to the long side of the bandage
- 2 Place the bandage 3 or 4 fingers below of the epicondyle (fig.A), aligning it with the extensor muscles of the hand
- 3 Insert the end part of the bandage into the appropriate ring; tighten the bandage and close it on itself (fig. B)
- 4 Check that the pressure pad is in the correct position
- 5 Should the bandage disturb too much when wearing it or cause a tingling sensation, slightly loosen the closure.

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Velcro® est une marque enregistrée de Velcro Industries B.V.