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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Medial epicondylitis strap providing compression

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medcal device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/ or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician. The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements. The product must be fitted with the utmost care to ensure its efficacy, tolerability and correct function. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. The manufacturer is no longer liable if the product is used or adapted inappropriately. The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used.

SELECTION/SIZE

Code	REF. 4320	
Size	S	L
Forearm circumference (cm)	< 27	> 27
Colour	beige	
fits right and left		

MAINTENANCE

⌘ Do not bleach ⌘ No chemical cleaning

⌘ Do not iron ⌘ Do not tumble-dry

🧼 Washing instructions: Remove the pads. Hand wash the bandage, the cover of the cooling pad and the pressure pad at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Severe and chronic epicondylitis and medial epicondylitis
- Muscle contractions
- Preventing inflammation during sports

CONTRAINDICATIONS

There are no known contraindications at this time.

FEATURES AND MATERIALS

- Double-lined neoprene with leatherette support
- Velcro® fastening with a nylon loop
- Non-slip fabric pressure pad

PUTTING ON THE APPLIANCE

- Attach the pressure pad on the bandage in the area Velcro® provided; the pressure pad has to be placed with the long part perpendicular to the long side of the bandage
- Place the bandage 3 or 4 fingers below of the epicondyle (fig.A), aligning it with the extensor muscles of the hand
- Insert the end part of the bandage into the appropriate ring; tighten the bandage and close it on itself (fig. B)
- Check that the pressure pad is in the correct position
- Should the bandage disturb too much when wearing it or cause a tingling sensation, slightly loosen the closure.

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Velcro® è un marchio di fabbrica registrato di Velcro Industries B.V.