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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM

REF.82300 Neoprene wrist wrap
REF.82600 Neoprene wrist and thumb wrap

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/ EEC and 2007/47/EC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. If in doubt about how to apply, contact an orthopaedic technician. Please read the composition of the product on the inside label carefully. You should not wear the device in the vicinity of an open flame. People prone to dermatitis or allergies to products resulting from the neoprene must not use this brace.

WARNINGS

The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician. The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements. The product must be fitted with the utmost care to ensure its efficacy, tolerability and correct function. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. The manufacturer is no longer liable if the product is used or adapted inappropriately. The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used.

SELECTION/SIZE

Code	REF.82300 / REF.82600
Size	One size
Color: beige · fits right and left	

MAINTENANCE

- ☒ Do not bleach
- ☒ No chemical cleaning
- ☒ Do not iron
- ☒ Do not tumble-dry
- ☒ Washing instructions: hand wash at max. 30° C with neutral soap (a sponge is recommended); do not dry in the vicinity of heat sources.

Do not dispose of the product or any of its components into the environment.

FEATURES AND MATERIALS

- Object 82300**
- Double-lined neoprene design
 - Velcro® fastening

- Object 82600**
- Double-lined neoprene design
 - Stabilising tubular thumb splint
 - Velcro® fastening

INDICATIONS

- Object 82300**
- Aftercare for wrist sprains
 - Tendonitis and tenosynovitis
 - Inflammatory degenerative diseases of the wrist

- Object 82600**
- Aftercare for a sprained thumb
 - Acute trapeziometacarpal arthritis
 - Arthrosis of the thumb

CONTRAINDICATIONS

Currently no known.

PUTTING ON THE APPLIANCE

- 1 Slide the thumb through the opening provided (A).
- 2 Hold the end of the wrist wrap firmly on the inside of the wrist (B).
- 3 Wrap the part with Velcro® around the wrist and fasten in place with the strap (C).

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