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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

**REF.82300 Neoprene wrist support**  
**REF.82600 Neoprene wrist support with thumb**

**DECLARATION OF CONFORMITY**

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

**TRADEMARKS OF MATERIALS**

Velcro® is a registered trademark of Velcro Industries B.V.

**SAFETY PRECAUTIONS**

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

**WARNINGS**

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

**SELECTION/SIZE**

Code	REF.82300 / REF.82600
Size	Universal
Color: beige · fits right and left	

**MAINTENANCE**

- Do not bleach
- Do not iron
- Washing instructions: hand wash at max. 30° C with neutral soap (a sponge is recommended); do not dry in the vicinity of heat sources.
- No chemical cleaning
- Do not tumble-dry

Do not dispose of the product or any of its components into the environment.

**INDICATIONS**

- Object 82300**
- Wrist sprains
  - Tendonitis and tenosynovitis
  - Inflammatory degenerative wrist diseases

**Object 82600**

- Sprained thumb
- Severe carpometacarpal osteoarthritis
- Thumb arthrosis

**CONTRAINDICATIONS**

Currently no known.

**FEATURES AND MATERIALS**

- Object 82300**
- Double-lined neoprene
  - Velcro® fastening system

**Object 82600**

- Double-lined neoprene
- Stay and stability of thumb
- Velcro® fastening system

**PUTTING ON THE APPLIANCE**

- 1 Slide the thumb through the opening provided (A).
- 2 Hold the end of the wrist wrap firmly on the inside of the wrist (B).
- 3 Wrap the part with Velcro® around the wrist and fasten in place with the strap (C).

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