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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Neoprene wrist brace with stay

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. People prone to dermatitis or allergies to products resulting from the neoprene must not use this brace. Do not apply in direct contact with open wounds.

WARNINGS




It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.8248			
Size	S	M	L	universal
Wrist circumf. cm	10/14	15/18	19/24	-
Length cm	18	23	23	18

Colour: black / specify right or left

MAINTENANCE

-  Do not bleach
-  No chemical cleaning
-  Do not iron
-  Do not tumble-dry
-  Washing instructions: Remove the palm splint from its seat. Hand wash in warm water with mild soap; rinse it thoroughly. Allow to dry away from heat sources

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Wrist and carpal sprains
- Tendinitis wrist
- Arthrosis and arthritis of wrist
- After wrist fractures
- Post-surgery wrist treatment

FEATURES AND MATERIALS

- Neoprene lined with antimicrobial nylon mesh material
- Flexible aluminium palm brace
- Velcro® fastening system

CONTRAINDICATIONS

Currently no known.

PUTTING ON THE APPLIANCE

(initial fitting only) If necessary, slide the metal stay out of its pocket (A) and shape it to fit the patient's wrist, bearing in mind that it should be positioned on the palm side of the hand.

- 1 Undo all the Velcro straps.
- 2 Rest the part of the brace containing the stay (A) on the palm of the hand with the stay pocket facing outwards.
- 3 Fasten the brace with the Velcro straps for the time being.
- 4 Secure the brace in position by first tightening the Velcro® straps on the forearm (B) and then the fastener for the thumb (C).

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