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37% Polyurethan (PU) - 29% Polyamid (PA) - 24% Polyester (PL) - 10% Sonstige Fasern (AF).

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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Long-size wrist immobilizer with a rigid and malleable stay

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the instructions of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of fire flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor, physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient.
To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately and, in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF. 90100		
Size	S	M	L
Wrist circ. cm	13/18	18/23	23/28
Length cm	21	21	21
Color	black		

Fits right and left

MAINTENANCE

- ☒ Do not bleach ☒ No chemical cleaning
☒ Do not iron ☒ Do not tumble-dry
☒ Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment

INDICATIONS

- Wrist distortion
- Compound fracture of the forearm of radius and ulna
- Wrist tendinitis
- Arthritis and arthrosis
- Consequences of fracture of the wrist or after early removal of the plaster

CONTRAINDICATIONS

There are no known contraindications at this time.

CHARACTERISTICS AND MATERIALS

- Fully welded, no seams
- Low profile
- Structure:
 - external: inextensible velveteen
 - padding: made of MTP next to the splints
 - internal (contact with skin): soft touch material (skin-compatible, comfortable, breathable, ladder-proof and antistatic)
- Thumb hole bordered in soft touch material for greater comfort
- Aluminum palm splint, preformed and removable on medical prescription
- Reinforcing fabric back of hand splints
- Velcro® fastening system
- Micro-velcro® hook to shorten the hand strap

PUTTING ON THE APPLIANCE

ONLY FOR THE FIRST USE (fig. A)

Pull the aluminum palm splint out of the orthesis and shape it to fit the patient's wrist. Reinsert the splint in the special pocket.

- 1 Open all the Velcro® straps (closing them again temporarily on themselves). If necessary open the orthesis fully and position it in such a way that the metal splint is located on the palm of the hand and the thumb is inserted in the appropriate opening (fig. B).
- 2 If it is not necessary to fully open the orthesis, put your hand in the wrist brace so that the aluminum palm splint is located on the palm of the hand and the thumb is inserted in the appropriate opening. (fig. C).
- 3 How to block the wrist brace: Tighten the central strap and fasten it using the Velcro® (fig. D).
- 4 Tighten the remaining straps and fasten them using the Velcro®: the hand strap (fig. E) and the forearm strap (fig. F).
- 5 If the hand strap is too long, it is possible to cut it along the special welding and fasten it using the back to back Velcro® hook included in the packaging.

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