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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Adjustable wrist-hand support brace

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.P1341 R · REF.P1341 L		
Size	S	M	L
Palm width (metacarpus) cm	5,5/7	7/8,5	8,5/10
Brace length cm	30	34	37
Color	pale blue		

specify right or left

MAINTENANCE

-  Do not bleach  No chemical cleaning
-  Do not iron  Do not tumble-dry
-  Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment

INDICATIONS

- Both conservative and surgical treatment of meta-distal epiphysis of forearm, fractures of carpus (except scaphoid), metacarpal fractures
- Effects of the fractures
- Hand and wrist tendinopathy
- Effects of flexors and finger extensors/tenorrhaphy
- Radio-carpal arthritis
- Neurologic deficits with subsequent paralysis

CONTRAINDICATIONS

- Mild and serious spasticity

CHARACTERISTICS AND MATERIALS

- Adjustable structure made of 2 metals with different stiffness
- Velcro® fastening system
- Ergonomic shape made of silver fiber in contact with the body

PUTTING ON THE APPLIANCE

DONNING THE BRACE FOR THE FIRST TIME WITH A DOCTOR TECHNICIAN

- Undo all the straps and fasten them back on themselves for the time being.
- Place the brace on the arm and hand with the palm facing upwards to check which alterations need to be made (fig. A).
- Adapt the brace to fit the patient and according to medical instructions by shaping the metal core:
 - shape the wrist area (this can be made easier by using a rounded surface for leverage) (fig. B)
 - adapt the forearm section (fig. C)
 - shape the area around the hand and fingers (fig. D).
- If necessary, the wrist strap can be shortened before fitting the brace on the patient by removing the Y-shaped Velcro® fastening, cutting the strap to the desired length, and repositioning the Y-shaped Velcro® fastening.
- Place the brace on the arm and hand with the palm facing upwards (fig. E) and for convenience first fasten the strap nearest the elbow by passing it through the ring provided and fastening it back on itself with the Velcro® (fig. F).
- Next fasten the remaining straps with Velcro®, namely the straps for the hand and fingers and the wrist strap (fig. G-H-I).

SUBSEQUENT WEAR

- Undo all the straps and fasten them back on themselves for the time being.
- Place the brace on the arm and hand with the palm facing upwards and for convenience first fasten the strap nearest the elbow by passing it through the ring provided and fastening it back on itself with the Velcro® (fig. F).
- Next fasten the remaining straps with Velcro®, namely the straps for the hand and fingers and the wrist strap (fig. G-H-I).

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