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REF. 5006

Thumb splint

DESIGNED AND DEVELOPED IN SWITZERLAND

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EN – The descriptions and images in this document are for illustrative and commercial purposes only. Orthoservice reserves the right to make modifications without prior notice based on their needs.

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Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13

info@orthoservice.de · www.orthoservice.de

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Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785

info@orthoservice.pl · www.orthoservice.pl

OSUK Ltd

Units 47&50 North East BIC

Wearfield

Sunderland Enterprise Park East

Sunderland SR5 2TA

Tel. +44 (0) 191 51 66 220

info@orthoservice.uk · www.orthoservice.uk

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
 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Thumb splint

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation.

In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.


SELECTION/SIZE

Code	REF. 5006			
Size	S	M	L	XL
Wrist circumference in cm	13/17	17/19	19/23	23/28
Colour	beige			

specify right or left

MAINTENANCE

 Do not bleach  No chemical cleaning
 Do not iron  Do not tumble-dry

 Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources. Plastic splint: Clean with a sponge soaked in handwarm water (30° C) and neutral soap. Dry with a cloth.

- Do not dispose of the product or any of its components into the environment
- Contact an orthopaedic technician for the replacement of used parts

INDICATIONS

- Rhizarthrosis
- Thumb tendinitis
- Tendinopathy of De Quervain
- Metacarpal/phalange sprain
- Steiner's lesion
- Effects of Bennet's fracture
- After trapezium plastic surgery and thumb surgery

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Smooth polypropylene structure
- Velcro® fastening system
- Equipped with antibacterial protective glove wear before brace
- Thermoformable

PUTTING ON THE APPLIANCE

- Put on the protective glove (A).
- Slip the plastic splint (B) on the thumb. If necessary, form it in order to adapt it to the thumb or choose another size.
- Tighten at the wrist with the Velcro® strap (C).

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