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
REF.2001

Finger splint (distal phalanx)

DESIGNED AND DEVELOPED IN SWITZERLAND

ORTHOSERVICE RO+TEN
Take care feel better

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EN – The descriptions and images in this document are for illustrative and commercial purposes only. Orthoservice reserves the right to make modifications without prior notice based on their needs.

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OSUK Ltd
Units 47&50 North East BIC
Wearfield
Sunderland Enterprise Park East
Sunderland SR5 2TA
Tel +44 (0) 191 51 66 220
info@orthoservice.uk · www.orthoservice.uk

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 0362 18 51 200 · Fax 0039 0362 18 58 065
info@roplusten.com · www.roplusten.com

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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Finger splint (distal phalanx)

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician.

Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames.

Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient.

To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation.

The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation.

In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.2001
Available sizes (For sizes, check the picture)	1 - 2 - 3 - 4 - 5 - 5 1/2 - 6 - 7
Colour	transparent

MAINTENANCE

- ⌘ Do not bleach
- ⌘ No chemical cleaning
- ⌘ Do not iron
- ⌘ Do not tumble-dry

Washing instructions: wipe the leather parts with a sponge dipped in lukewarm water and mild soap. Leave to dry away from direct heat.

Do not dispose of the product or any of its components into the environment

INDICATIONS

- Injury to extensor tendon
- Segond's lesion
- Fractures of the distal phalanx
- Sprain of distal interphalangeal joint
- Bruised nails
- Post-surgery of extensor tendon

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Clear plastic
- Rigid plastic and transpiration holes

PUTTING ON THE APPLIANCE

1 Slip the finger into the protector, making sure, the nail coincides with the free upper part.

2 Use tape or velcro to fasten the protector to the finger.

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