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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

REF.H1133 Thigh support in Aer-OS fabric

REF.A1134 Calf support in Aer-OS fabric

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician. The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements. The product must be fitted with the utmost care to ensure its efficacy, tolerability and correct function. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. The manufacturer is no longer liable if the product is used or adapted inappropriately. The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used.

SELECTION/SIZE

Item	REF. H1133				
Size	S	M	L	XL	XXL
Circ. middle of thigh cm	48/52	52/56	56/61	61/66	66/72
Colour	black				
fits right and left					

Item	REF. A1134			
Size	S	M	L	XL
Circ. calf cm	33/35	35/38	38/43	43/48
Colour	black			
fits right and left				

MAINTENANCE

-  Do not bleach
-  Do not iron
-  Do not tumble-dry
-  Washing instructions: Hand wash in lukewarm water (30°C) with mild soap; rinse thoroughly. Leave to dry away from direct heat.

Dispose of this product and its parts responsibly

INDICATIONS

MioFIT33

- Muscular contractions and tears of the thigh
- Aid to resuming sporting activity after muscle trauma
- Prevention of bruising

MioFIT34

- Muscular contractions and tears of the calf
- Aid to resuming sporting activity after muscle trauma

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

MioFIT33 – MioFIT34

Made from breathable Aer-OS fabric

PUTTING ON THE APPLIANCE

MioFIT33

- 1 Pull the support over the thigh to below the inner pelvic area; the orange label must face upwards (fig. A).
- NOTE:** The orange profile must be on the front

MioFIT34

- 1 Pull the support over the leg to below the knee cap so that the calf area is completely protected (as shown in the figure); the orange label must face upwards (fig. B).
- NOTE:** The orange profile must be on the back.

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