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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Neoprene thigh support

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS



It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.8233				
Size	S	M	L	XL	XXL
Circumference at mid-thigh cm	48/51	52/55	56/60	61/65	66/71
Colour	blue				

Fits right and left

MAINTENANCE

-  Do not bleach
-  No chemical cleaning
-  Do not iron
-  Do not tumble-dry
-  Washing instructions: hand wash at max. 30° C with neutral soap (a sponge is recommended); do not dry in the vicinity of heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Muscular contractions and tears of the thigh
- Aid to resuming sporting activity after trauma
- Prevention of bruising

CONTRAINDICATIONS

Currently no known.

FEATURES AND MATERIALS

- Double-lined orthopedic neoprene
- Strengthened neoprene patch
- Washable

PUTTING ON THE APPLIANCE

- 1 Insert thigh support up to the inner groin area.
- 2 Place the patch on the front side.

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