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Innovative tensioning system with pulley technology

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EASY PULLING

BREATHABLE

PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Semi-rigid corset for the sacroiliac joint

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The response to the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.T1032 R - T1032 L
Size	Regular Long
Hip circumference cm	< 89 89/120
Colour	black

MAINTENANCE

☒ Do not bleach ☒ No chemical cleaning
☒ Do not iron ☒ Do not tumble-dry
☒ Washing instructions: Hand wash at max. 30°C. Hand wash in lukewarm water with mild soap (at max. 30°C), rinse the device or any of its components. Dry away from heat sources.

Do not dispose of the device or any of its components.

SCOPE OF DELIVERY

The Symplus includes:

- 1 x Belt (long REF. T1032 R or short version REF. T1032 L) (fig. 1).
- 2 x padded pelottes with Velcro® closure (fig. 2).

INDICATIONS

- Sacroiliac joint (SI) pain or dysfunction
- Arthrosis of the sacroiliac joint
- Symphysis Pubis Dysfunction (SPD) (after childbirth). Note: in this case do not use the padded pads

CONTRAINDICATIONS

- In case of Symphysis Pubis Dysfunction (SPD) (after childbirth) do not use the padded pelottes

CHARACTERISTICS AND MATERIALS

- Innovative tensioning system with pulley technology for optimum stabilisation and for simple application minimal effort
- firm, non-deformable and breathable fabric for safe hold, with edges in anti-irritation material. Particularly light
- The large Velcro® fastening enables the tensioning strap to be secured in a simple way and to the desired compression.
- Reinforced paravertebral area
- Two padded pads that can be positioned differently thanks to the surfaces of Velcro® closure, designed to allow targeted compression

PUTTING ON THE APPLIANCE

ADAPTING FOR THE DOCTOR / TECHNICIAN

- If necessary, position the padded pelottes, as prescribed, using the Velcro® closure on the belt (fig. A).
- Open the belt and provisionally attach the Velcro® closure element to the side of the belt.
- Position the belt in such a way that the labels „right/left” are on the correct sides and are pointing upwards.
- Put on the belt and secure using the central Velcro® closure (fig. B).
- Optimise the tension of the corsets using the end piece (fig. C) and secure via the Velcro® closure on the appropriate Velcro® surface (fig. D).
- Make sure that the reinforced section at the back are positioned centrally/symmetrically to the spine.

PUTTING ON THE APPLIANCE FOR THE PATIENT

- Open the belt and provisionally attach the Velcro® closure element to the side of the belt.
- Position the belt in such a way that the labels „right/left” are on the correct sides and are pointing upwards.
- Put on the belt and secure using the central Velcro® closure (fig. B).
- Optimise the tension of the corsets using the end piece (fig. C) and secure via the Velcro® closure on the appropriate Velcro® surface (fig. D).
- Make sure that the reinforced section at the back are positioned centrally/symmetrically to the spine.

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