LÖZÜ, UZMANIĞI İZLASIET UN SAGLABÄJİET INSTRUKCIJAS

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  **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**
Semi-rigid corset for the sacroiliac joint

DECLARATION OF CONFORMITY
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.
TRADEMARKS OF MATERIALS
Velcro® is a registered trademark of Velcro Industries B.V.
SAFETY PRECAUTIONS
We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not too tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

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MAINTENANCE
⚠ Do not bleach ⚠ No chemical cleaning
⚠ Do not iron ⚠ Do not tumble-dry
Washing instructions: Hand wash at max. 30°C. Hand wash in lukewarm water with mild soap (at max. 30°C); rinse thoroughly. Dry away from heat sources.
Do not dispose of the device or any of its components.

SCOPE OF DELIVERY
The Symplus includes:
• 1 x Belt (long REF. T1032 R or short version REF. T1032 L) (fig. 1).
• 2 x padded pelottes with Velcro® closure (fig. 2).

INDICATIONS
• Sacroiliac joint (SI) pain or dysfunction
• Arthrosis of the sacroiliac joint
• Symphysis Pubis Dysfunction (SPD) (after childbirth). Note: in this case do not use the padded pads

CONTRAINDICATIONS
• In case of Symphysis Pubis Dysfunction (SPD) (after childbirth) do not use the padded pelottes

CHARACTERISTICS AND MATERIALS
• Innovative tensing system with pulley technology for optimum stabilisation and for simple application minimal effort
• Firm, non-deformable and breathable fabric for safe hold, with edges in anti-irritation material. Particularly light
• The large Velcro® fastening enables the tensing strap to be secured in a simple way and to the desired compression.
• Reinforced paravertebral area
• Two padded pads that can be positioned differently thanks to the surfaces of Velcro® closure, designed to allow targeted compression

PUTTING ON THE APPLIANCE
ADAPTING FOR THE DOCTOR / TECHNICIAN
1. If necessary, position the padded pelottes, as prescribed, using the Velcro® closure on the belt (fig. A).
2. Open the belt and provisionally attach the Velcro® closure element to the side of the belt.
3. Position the belt in such a way that the labels „right/left“ are on the correct sides and are pointing upwards.
4. Open the belt and secure using the central Velcro® closure (fig. B).
5. Optimise the tension of the corsets using the end piece (fig. C) and secure via the Velcro® closure on the appropriate Velcro® surface (fig. D).
6. Make sure that the reinforced section at the back are positioned centrally/symmetrically to the spine.

PUTTING ON THE APPLIANCE FOR THE PATIENT
1. Open the belt and provisionally attach the Velcro® closure element to the side of the belt.
2. Position the belt in such a way that the labels „right/left“ are on the correct sides and are pointing upwards.
3. Put on the belt and secure using the central Velcro® closure (fig. B).
4. Optimise the tension of the corsets using the end piece (fig. C) and secure via the Velcro® closure on the appropriate Velcro® surface (fig. D).
5. Make sure that the reinforced section at the back are positioned centrally/symmetrically to the spine.

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